

Europäisch s Patentamt

**European Patent Office** 

Office européen des brevets



(11) EP 0 668 914 B1

(12)

# **EUROPEAN PATENT SPECIFICATION**

- (45) Date of publication and mention of the grant of the patent: 16.08.2000 Bulletin 2000/33
- (21) Application number: 93915279.9
- (22) Date of filing: 09.06.1993

- (51) Int. Cl.<sup>7</sup>: **C12N 15/27**, C12P 21/02, C07K 14/53, C12N 15/12, G01N 33/68
- (86) International application number: PCT/US93/05548
- (87) International publication number:
   WO 93/25687 (23.12.1993 Gazette 1993/30)

# (54) CRYSTALLIZATION OF M-CSF

KRISTALLISIERUNG VON M-CSF CRISTALLISATION DU M-CSF

- (84) Designated Contracting States:

  AT BE CH DE DK ES FR GB GR IE IT LI LU MC NL
  PT SE
- (30) Priority: 09.06.1992 US 896512
- (43) Date of publication of application: 30.08.1995 Bulletin 1995/35
- (60) Divisional application: 99108304.9 / 0 955 365
- (73) Proprietors:
  - CHIRON CORPORATION
     Emeryville California 94608-2916 (US)
  - The Regents of the University of California Oakland, CA 94607-5200 (US)
- (72) Inventors:
  - PANDIT, Jayvardhan Mystic, CT 06355 (US)
  - JANCARIK, Jarmila Walnut Creek, CA 94596 (US)
  - KIM, Sung-Hou Moraga, CA 94556 (US)
  - KOTHS, Kirston El Cerrito, CA 94530 (US)
  - HALENBECK, Robert San Rafael, CA 94901 (US)
  - FEAR, Anna, Lisa
     Oakland, CA 94609 (US)
  - TAYLOR, Eric Oakland, CA 94609 (US)

- YAMAMOTO, Ralph Martinez, CA 94553 (US)
- BOHM, Andrew Berkeley, CA 94704 (US)
- (74) Representative:
  Ritter, Stephen David et al
  Mathys & Squire
  100 Grays Inn Road
  London WC1X 8AL (GB)
- (56) References cited:

EP-A- 0 344 796 WO-A-88/08003

EP-A- 0 402 070 WO-A-92/21029

- R. K. SCOPES 'protein purification advanced texts in chemistry' 1987, SECOND EDITION. PAGES 297-321. SPRINGER-VERLAG, NEW YORK, US.
- JOURNAL OF BIOLOGICAL CHEMISTRY. vol. 251, no. 20, 25 October 1976, BALTIMORE US pages 6300 - 6303 A. MCPHERSON, JR 'Crystallization of proteins from polyethylene glycol' cited in the application
- SCIENCE vol. 258, 20 November 1992, LANCASTER, PA US pages 1358 - 1362 J. PANDIT ET AL 'Three-dimensional structure of dimeric human recombinant Macrophage Colony-Stimulating Factor'
- SPORN M.B. ET ROBERTS A.B. 'peptide growth factors and their receptors, handbook of experimental pharmacology vol. 1; vol. 95/1. pages 667-698.' 1990, SPINGER VERLAG, DE, BERLIN.

Note: Within nine months from the publication of the mention of the grant of the European patent, any person may give notice to the European Patent Office of opposition to the European patent granted. Notice of opposition shall be filed in a written reasoned statement. It shall not be deemed to have been filed until the opposition fee has been paid. (Art. 99(1) European Patent Convention).

- PROCEEDINGS OF THE NATIONAL ACADEMY
  OF SCIENCES OF USA. vol. 84, March 1987,
  WASHINGTON US pages 1268 1271 Y. G.
  YEUNG ET AL 'Purification of the colony
  stimulating factor 1 receptor and
  demonstration of its tyrosin kinase activity'
- JOURNAL OF MOLECULAR BIOLOGY vol. 205, no. 4, 1989, pages 783 - 785 J. M. LA LONDE ET AL 'Crystallization and preliminary X-ray studies of recombinant human granulocyte - macrophage colony stimulating factor'
- SCIENCE vol. 254, 8 November 1991, LANCASTER, PA US page 821825 B. C.
   CUNNINGHAM ET AL 'Dimerization of the extracellular domain of the human growth hormone receptor by a single hormone molecule'
- JONES, N.D. ET AL.: "Crystallization of authentic recombinant human growth hormone", BIOTECHNOLOGY, vol. 5, pages 499 to 500, May 1987

#### Description

#### **CRYSTALLIZATION OF M-CSF**

[0001] Work described herein was funded with the United States Government support. The United States Government has certain rights in inventions arising as part of that work. This application is a continuation-in-part of U.S. Serial No. 07/896,512 filed June 9, 1992.

#### FIELD OF THE INVENTION

10

[0002] The present invention relates in general to crystalline compositions of macrophage colony stimulating factor "M-CSF" and in particular to methods for the use of structural information (including X-ray diffraction patterns) of crystalline M-CSF for agonist and antagonist production, as well as assays for detection of same.

#### BACKGROUND OF THE INVENTION

[0003] Monocyte-macrophage colony-stimulating factor is produced by a variety of cells, including macrophages, endothelial cells and fibroblasts (see, Ralph et al., "The Molecular and Biological Properties of the Human and Murine Members of the CSF-1 Family" in Molecular Basis of Lymphokine Action, Humana Press, Inc., (1987),). M-CSF is composed of two "monomer" polypeptides, which form a biologically active dimeric M-CSF protein (hereinafter referred to as "M-CSF dimer"). M-CSF belongs to a group of biological agonists that promote the production of blood cells. Specifically, it acts as a growth and differentiation factor for bone marrow progenitor cells of the mononuclear phagocyte lineage. Further, M-CSF stimulates the proliferation and function of mature macrophages via specific receptors on responding cells. In clinical trials M-CSF has shown promise as a pharmaceutical agent in the correction of blood cell deficiencies arising as a side-effect of chemotherapy or radiation therapy for cancer and may be beneficial in treating fungal infections associated with bone marrow transplants. M-CSF may also play significant biological roles in pregnancy, uveitis, and atherosclerosis. Development of M-CSF agonists or antagonists may prove to be of value in modifying the biological events involved in these conditions.

[0004] M-CSF exists in at least three mature forms: short (M-CSFα), intermediate (M-CSFγ), and long (M-CSFβ). Mature M-CSF is defined as including polypeptide sequences contained within secreted M-CSF following amino terminus processing to remove leader sequences and carboxyl terminus processing to remove domains including a putative transmembrane region. The variations in the three mature forms are due to alternative mRNA splicing (see, Cerretti et al. Molecular Immunology, 25:761 (1988)). The three forms of M-CSF are translated from different mRNA precursors, which encode polypeptide monomers of 256 to 554 amino acids, having a 32 amino acid signal sequence at the amino terminal and a putative transmembrane region of approximately 23 amino acids near the carboxyl terminal. The precursor peptides are subsequently processed by amino terminal and carboxyl terminal proteolytic cleavages to release mature M-CSF. Residues 1-149 of all three mature forms of M-CSF are identical and are believed to contain sequences essential for biological activity of M-CSF. In vivo M-CSF monomers are dimerized via disulfide-linkage and are glycosylated. Some, if not all, forms of M-CSF can be recovered in membrane-associated form. Such membrane-bound M-CSF may be cleaved to release secreted M-CSF. Membrane associated M-CSF is believed to have biological activity similar to M-CSF, but may have other activities including cell-cell association or activation.

[0005] Polypeptides, including the M-CSFs, have a three-dimensional structure determined by the primary amino acid sequence and the environment surrounding the polypeptide. This three-dimensional structure establishes the polypeptide's activity, stability, binding affinity, binding specificity, and other biochemical attributes. Thus, a knowledge of a protein's three-dimensional structure can provide much guidance in designing agents that mimic, inhibit, or improve its biological activity in soluble or membrane bound forms.

[0006] The three-dimensional structure of a polypeptide may be determined in a number of ways. Many of the most precise methods employ X-ray crystallography (for a general review, see, Van Holde, *Physical Biochemistry, Prentice-Hall, N.J.* pp. 221-239, (1971), which is incorporated herein by reference). This technique relies on the ability of crystalline lattices to diffract X-rays or other forms of radiation. Diffraction experiments suitable for determining the three-dimensional structure of macromolecules typically require high-quality crystals. Unfortunately, such crystals have been unavailable for M-CSF as well as many other proteins of interest. Thus, high-quality, diffracting crystals of M-CSF would assist the determination of its three-dimensional structure.

[0007] Various methods for preparing crystalline proteins and polypeptides are known in the art (see, for example, McPherson, et al. "Preparation and Analysis of Protein Crystals", A. McPherson, Robert E. Krieger Publishing Company, Malabar, Florida (1989); Weber, Advances in Protein Chemistry 41:1-36 (1991); U.S. Patent No. 4,833,233;). Although there are multiple approaches to crystallizing polypeptides, no single set of conditions provides a reasonable expectation of success, especially when the crystals must be suitable for X-ray dif-

fraction studies. Thus, in spite of significant research, many proteins remain uncrystallized.

[0008] In addition to providing structural information, crystalline polypeptides provide other advantages. For example, the crystallization process itself further purifies the polypeptide, and satisfies one of the classical criteria for homogeneity. In fact, crystallization frequently provides unparalleled purification quality, removing impurities that ar not removed by other purification methods such as HPLC, dialysis, conventional column chromatography, etc. Moreover, crystalline polypeptides are often stable at ambient temperatures and free of protease contamination and other degradation associated with solution storage. Crystalline polypeptides may also be useful as pharmaceutical preparations. Finally, crystallization techniques in general are largely free of problems such as denaturation associated with other stabilization methods (e.g. lyophilization). Thus, there exists a significant need for preparing M-CSF compositions in crystalline form and determining their three-dimensional structure. The present invention fulfills this and other needs. Once crystallization has been accomplished, crystallographic data provides useful structural information which may assist the design of peptides that may serve as agonists or antagonists. In addition, the crystal structure provides information useful to map, the receptor binding domain which could then be mimicked by a small non-peptide molecule which may serve as an antagonist or agonist.

# SUMMARY OF THE INVENTION

15

[0009] The present invention provides crystalline forms of M-CSF dimers. Preferably, the dimers are formed from polypeptides containing between 146 to 162 amino acids residues at or near the N-terminus of mature M-CSF (e.g.  $glu_1$   $glu_2$   $val_3$ ...). In a specific embodiment, the polypeptide includes residues 4 to 158 of mature M-CSF $\alpha$  polypeptide, preferably in the non-glycosylated form.

[0010] Another aspect of the invention provides a method of crystallizing an M-CSF. A preferred crystallization method according to the present invention includes the following steps: mixing a preselected, substantially pure M-CSF dimer and a precipitant to form an M-CSF mixture; precipitating crystalline M-CSF from the mixture; and isolating the M-CSF in crystalline form. In some specific embodiments, the precipitant contains polyethylene glycol. Other components such as ammonium sulfate and/or other ionic compounds may be added to the solution. It has been found by x-ray crystallography that M-CSF produced by the method of the present invention can crystallize into the P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub> space groups for example.

Variations of the crystallization method are also provided. For example, the step of precipitating crystals [0011] from the M-CSF mixture may involve equilibrating the M-CSF mixture with a second mixture. The second mixture is typically a solution that consists of a higher concentration of precipitant than the first M-CSF mixture. The step of equilibrating preferably consists of applying the M-CSF mixture to a surface and allowing the applied M-CSF mixture to come into equilibrium with a reservoir of the second mixture. In other embodiments, the step of precipitating M-CSF crystals is initiated by seeding the M-CSF mixture with seed crystals or altering the temperature of the M-CSF mixture. Another aspect of the invention involves identifying compounds that have structures that mimic a receptor binding region of the three-dimensional structure of M-CSF to varying degrees and can in many instances function as M-CSF agonists or antagonists. Compounds that interact with the receptor-binding region of M-CSF may be antagonists. The three-dimensional alpha-carbon coordinates of a truncated M-CSF dimer is presented in Appendix 1. In one embodiment of the present invention, the three-dimensional structure of M-CSF is obtained by first crystallizing an M-CSF dimer (having M-CSF receptor-binding residues) to form at least one M-CSF crystal. Next, a source of radiation is used for irradiating an M-CSF crystal to obtain a diffraction pattern of the M-CSF crystal. Finally, a three-dimensional structure of M-CSF is obtained from the diffraction pattern. In most embodiments, the three-dimensional structure includes an M-CSF receptor-binding region.

Structural features identified in accordance with the invention may be used for selecting candidate amino [0012]acid substitutions in a protein, based on structural information, and more particularly M-CSF, comprising determining the three-dimensional structure of M-CSF by the methods of the present invention; followed by determining the solvent accessible amino acid residues of the protein, determining which residues are not involved in dimer formation. Applying these criteria, amino acids in M-CSF which are solvent accessible and which are not involved in dimer formation are selected for substitution with non-conservative amino acids. Since M-CSFB has intrachain disulfide bonds involving cysteines 157 and/or 159, we believe the C-terminal region of M-CSF extends from the "rear" of the structure we have solv d, providing a variable-length "tether" for membrane-bound forms of M-CSF. Thus, the "front" or receptor-binding region of M-CSF is on the opposite side of the molecules, consisting of solvent-accessible residues in or near helices A, C, and D, including residues from about 6 to 26, 71 to 90, and 110 to 130, respectively, of native M-CSF. Preferred amino acids for substitution and preferred substituting amino acids include but are not limited to: H15  $\rightarrow$  A or L; Q79  $\rightarrow$ A or D; R86  $\rightarrow$  E or D; E115  $\rightarrow$  A; E41  $\rightarrow$  K or R; K93  $\rightarrow$  A or E; D99  $\rightarrow$  K or R; L55  $\rightarrow$  Q or N; S18  $\rightarrow$  A or K; Q20  $\rightarrow$ A or D; 175  $\rightarrow$  K or E; V78  $\rightarrow$  K or R; L85  $\rightarrow$  E or N; D69  $\rightarrow$  K or R; N70  $\rightarrow$  A or E; H9  $\rightarrow$  A or D; N63  $\rightarrow$  K or R; and  $T34 \rightarrow Q$  or K. Most preferred are those substitutions give rise to novel M-CSF agonists and M-CSF antagonists. Additionally, the present invention is also directed to a method for producing antagonists and agonists by substituting at least

one and preferably fewer than 5 solvent accessible residues per M-CSF monomer.

[0013] The invention also enables the production of heterodimeric M-CSF in which only one subunit contains substituted solvent accessible amino acids involved in signal transduction and to heterodimeric M-CSF in which each subunit contains different substituted solvent accessible amino acids involved in signal transduction. The present invention is also enables the production of M-CSF having amino acid substitutions which do not impair binding to the M-CSF receptor. Screening for agonists and antagonists is then accomplished using bioassays and receptor binding assays using methods well known in the art, including those described in the Examples below.

[0014] In addition the invention is directed to an isolated, purified, soluble and functional M-CSF receptor. The present invention is also directed to a method for screening M-CSF agonists and antagonists using a soluble M-CSF receptor.

[0015] A further understanding of the present invention can be obtained by reference to the drawings and discussion of specific embodiments.

# **BRIEF DESCRIPTION OF THE FIGURES**

#### [0016]

15

30

35

40

Figure 1 is a section of a diffraction pattern of an M-CSF $\alpha$  crystal prepared according to the present invention; Figure 2 is a topology diagram showing the disulfide bonds in truncated dimeric M-CSF;

Figure 3 is a stereodiagram of the C-alpha backbone with every tenth residue labelled and with the non-crystallographic symmetry axis indicated by a dotted line;

Figures 4A and 4B present two views of a ribbon diagram highlighting the secondary structural elements of M-CSF. The cysteine residues have been represented by a ball-and stick model and the non-crystallographic symmetry axis is indicated by a dotted line;

Figures 5A-C illustrate size-exclusion HPLC analysis of the NΔ3CΔ158 M-CSF short clone homodimer (158) (Figure 5A), NΔ3CΔ221 C157S, C159S long clone homodimer (221F) (Figure 5B), and the short clone/long clone heterodimer (158/221F) (Figure 5C) and their corresponding biological activities;

Figures 6A-6C illustrates size exclusion HPLC (Figure 6A) and both non-reduced (Figure 6B) and reduced (Figure 6C) SDS-PAGE analysis of the preparative purification of M-CSF, more particularly, Figure 6A graphically illustrates the separation on Phenyl-HPLC size exclusion chromatography of the three species of M-CSF dimers of Figures 5A-5C, i.e., the 158 homodimer, the 221F homodimer and the 158/221 heterodimer, and indicates that the three absorbance peaks at 280 nm (solid line) correlate with M-CSF activity in U/ml x 10<sup>-6</sup> (dotted line); Figure 6B illustrates an SDS-PAGE analysis under non-reducing conditions of preparative purification of the 158/221F heterodimer (intermediate molecular weight species) relative to the 158 homodimer (lower molecular weight species) and the 221F homodimer (highest molecular weight species); Figure 6C illustrates an SDS-PAGE analysis under reducing conditions of the preparative purification of the 158/221F heterodimer (middle lanes) relative to the 158 homodimer (left lanes) and the 221F homodimer (right lanes): and

Figure 7 illustrates the competitive binding of M-CSF and M-CSF muteins to NFS60 cell M-CSF receptors. In Figure 7, competitive binding curves are shown for M-CSF $\alpha$  N $\Delta$ 3C $\Delta$ 158 (closed circles); M-CSF $\alpha$  N $\Delta$ 3C $\Delta$ 158 H9A, H15A/M-CSF $\beta$  N $\Delta$ 3C $\Delta$ 221 C157S, C159S heterodimer (closed squares); dimeric Q20A, V78KF mutein (open circles); and dimeric H9A, H15A mutein (open squares).

# DESCRIPTION OF THE PREFERRED EMBODIMENTS DEFINITIONS

5 [0017] As used herein "M-CSF polypeptide" refers to a human polypeptide having substantially the same amino acid sequence as the mature human M-CSFα, M-CSFβ, or M-CSFγ polypeptides described in Kawasaki et al, Science 230:291 (1985), Cerretti et al, Molecular Immunology, 25:761 (1988), or Ladner et al., EMBO Journal 6:2693 (1987). Such terminology reflects the understanding that the three mature M-CSFs have different amino acid sequences, as described above.

Certain modifications to the primary sequence of M-CSF can be made by deletion, addition, or alteration of the amino acids encoded by the DNA sequence without destroying the desired structure (e.g., the receptor binding ability of M-CSF) in accordance with well-known recombinant DNA techniques. Further, a skilled artisan will appreciate that individual amino acids may be substituted or modified by oxidation, reduction or other derivitization, and the polypeptide may be cleaved to obtain fragments that retain the active binding site and structural information. Such substitutions and alterations result in polypeptides having an amino acid sequence which falls within the definition of polypeptide "having substantially the same amino acid sequence as the mature M-CSFα, M-CSFβ, and M-CSFγ polypeptides."

[0019] For purposes of crystallization, preferred lengths of the M-CSFα, β or γ monomers are between about 145 and 180 amino acids (counting from the mature amino terminus), and more preferably between about 145 and 162

amino acids long. A specific monomer that may be present in a crystallizable dimer is M-CSF $\alpha$  and is N $\Delta$ 3M-CSF $\alpha$ C $\Delta$ 158 (3 amino acids are deleted from the amino terminus and the total length is 155 amino acids). All lengths are inclusive. As used herein the term "M-CSF $\alpha$  (4-158)" denotes an M-CSF having amino acid residues 4 to 158 of the mature, processed M-CSF $\alpha$  polypeptide. Other nomenclature designations for C-terminal and N-terminal truncations of native M-CSF are set forth in U.S. Patent No. 4,929,700.

[0020] Crystallizable glycosylation variants of the M-CSF polypeptides are included within the scope of this invention. These variants include polypeptides completely lacking in glycosylation and variants having at least one fewer glycosylated site than the mature forms, as well as variants in which the glycosylation pattern has been changed from the native forms. Also included are deglycosylated and unglycosylated amino acid sequence variants, as well as deglycosylated and unglycosylated M-CSF subunit having the mature amino acid sequence (see, U.S. Patent No. 5,032,626).

[0021] "M-CSF" dimer refers to two M-CSF polypeptide monomers that have dimerized M-CSF dimers may include two identical polypeptide monomers (homodimers) or two different polypeptide monomers (heterodimers such as an M-CSFα-M-CSFβ dimer, an M-CSF long clone and short clone dimer). M-CSF monomers may be convened to M-CSF dimers *in vitro* as described in U.S. Patent No. 4,929,700. Recombinantly expressed M-CSFs may also be variably glycosylated as they exist *in vivo*, partially glycosylated, or completely lacking in glycosylation (unglycosylated). Glycosylated M-CSFs may be produced *in vivo* with carbohydrate chains which may later be enzymatically deglycosylated *in vitro*.

[0022] Biologically active M-CSF exhibits a spectrum of activity understood in the art. For instance, M-CSF stimulates the proliferation and function of mature macrophages via specific receptors on responding cells. Further, M-CSF acts as a mononuclear phagocyte progemtor growth factor. The standard *in vitro* colony stimulating assay of Metcalf, *J. Cell Physiol.* 76:89 (1970) results primarily in the formation of macrophage colonies when M-CSF is applied to stem cells. Other biological assays are based on M-CSF induced proliferation of M-CSF dependent cells such as the NFS-60 cell line. As used herein "M-CSF having biological activity" refers to M-CSF, including fragments and sequence variants thereof as described herein; that exhibit an an-recognized spectrum of activity with respect to biological systems. Such M-CSF having biological activity will typically have certain structural attributes in common with those of the mature M-CSF forms such as receptor binding site tertiary structure.

[0023] Agonists are substances that exhibit greater activity *per se* than the native ligand while antagonists are substances that suppress, inhibit, or interfere with the biological activity of a native ligand. Agonists and antagonists may be produced by the methods of the present invention for use in the treatment of diseases in which M-CSF has been implicated either as a potential treatment (*e.g.*, for treating blood cell deficiencies arising as a side effect of chemotherapy treating fungal infection associated with bone marrow transplants and others) or as having a role in the pathogenesis of the disease (*e.g.*, ovarian cancer, uveitis, atherosclerosis).

[0024] Crystallization of M-CSF species in accordance with the present invention includes four general steps: expression, purification, crystallization and isolation.

#### Expression of Recombinant M-CSF

35

45

[0025] M-CSF crystallization requires an abundant source of M-CSF that may be isolated in a relatively homogeneous form. A variety of expression systems and hosts are suitable for the expression of M-CSF and will be readily apparent to one of skill in the art. Because of the variability of glycosylation and other post-transnational modifications present in M-CSF produced in certain eukaryotic hosts, expression in *E. coli* may provide M-CSF with advantageous properties with regard to crystallization. Typical *in vitro* M-CSF expression systems are described in U.S. Patent No. 4,929,700, for example.

[0026] For use in the present invention, a variety of M-CSF polypeptides can also be readily designed and manufactured utilizing recombinant DNA techniques well known to those skilled in the art. For example, the M-CSF amino acid sequence can vary from the naturally occurring sequence at the primary structure level by amino acid substitutions, insertions, deletions, and the like. These modifications can be used in a number of combinations to produce the final modified polypeptide chain. The present invention is useful for crystallizing such polypeptides and dimers thereof.

[0027] In general, modifications of the genes encoding the M-CSF polypeptide are readily accomplished by a variety of well-known techniques, such as site-directed mutagenesis (see, Gillman and Smith, Gene 8:81-97 (1979) and Roberts, S. et al., Nature 328:731-734 (1987) and U.S. Patent No. 5,032,676,). Most modifications are evaluated by screening in a suitable assay for the desired characteristic. For instance, a change in the M-CSF receptor-binding character of the polypeptide can be detected by competitive assays with an appropriate reference polypeptides or by the bioassays described in U.S. Patent No. 4,847,201, which is incorporated herein by reference.

[0028] Insertional variants of the present invention are those in which one or more amino acid residues are introduced into a predetermined site in the M-CSF. For instance, insertional variants can be fusions of heterologous proteins or polypeptides to the amino or carboxyl terminus of the subunits. Substitutional variants are those in which at least one residue has been removed and a different residue inserted in its place. Nonnatural amino acids (i.e., amino acids not

normally found in native proteins), as well as isosteric analogs (amino acid or otherwise) are also suitable for use in this invention. Examples of suitable substitutions are well known in the art, such as the Glu->Asp, Ser->Cys, and Cys->Ser, His-> alanine for example. Another class of variants are deletional variants, which are characterized by the removal of one or more amino acid residues from the M-CSF.

- Other variants of the present invention may be produced by chemically modifying amino acids of the native protein (e.g., diethylpyrocarbonate treatment which modifies histidine residues). Preferred or chemical modifications which are specific for certain amino acid side chains. Specificity may also be achieved by blocking other side chains with antibodies directed to the side chains to be protected. Chemical modification includes such reactions as oxidation, reduction, amidation, deamidation, or substitution of bulky groups such as polysaccharides or polyethylene glycol (see e.g., U.S. Patent No. 4,179.337 and WO91/21029).
  - [0030] Exemplary modifications include the modification of lysinyl and amino terminal residues by reaction with succinic or other carboxylic acid anhydrides. Modification with these agents has the effect of reversing the charge of the lysinyl residues. Other suitable reagents for modifying amino-containing residues include imidoesters such as methyl picolinimidate; pyridoxal phosphate; pyridoxal chloroborohydride; trinitrobenzenesulfonic acid; O-methylisourea, 2,4-pentanedione; and transaminase-catalyzed reaction with glyoxylate, and N-hydroxysuccinamide esters of polyethylenene glycol or other bulky substitutions.
  - [0031] Arginyl residues may be modified by reaction with a number of reagents, including phenylglyoxal, 2,3-butanedione, 1,2-cyclohexanedione, and ninhydrin. Modification of arginine residues requires that the reaction be performed in alkaline conditions because of the high  $pK_a$  of the guanidine functional group. Furthermore, these reagents may react with the groups of lysine as well as the arginine epsilon-amino group.
  - [0032] Tyrosyl residues may also be modified with particular interest in introducing spectral labels into tyrosyl residues by reaction with aromatic diazonium compounds or tetranitromethane are used to form O-acetyl tyrosyl species and 3-nitro derivatives, respectively. Tyrosyl residues may also be iodinated using <sup>125</sup>I or <sup>131</sup>I to prepare labeled proteins for use in radioimmunoassay.
- 25 [0033] Carboxyl side groups (aspartyl or glutamyl) may be selectively modified by reaction with carbodiimides (R-N=C=N-R¹), where R and R¹ are different alkyl groups, such as 1-cyclohexyl-3-(2-morpholinyl-4-ethyl) carbodiimide or 1-ethyl-3-(4-azonia-4,4-dimethylpentyl) carbodiimide. Furthermore, aspartyl and glutamyl residues are convened to asparaginyl and glutaminyl residues by reaction with ammonium ions.
- [0034] Conversely, glutaminyl and asparaginyl residues may be deamidated to the corresponding glutamyl and aspartyl residues, respectively, under mildly acidic conditions. Either form of these residues falls within the scope of this invention.
  - [0035] Other modifications include hydroxylation of proline and lysine, phosphorylation of hydroxyl groups of seryl or threonyl residues, methylation of the  $\alpha$ -amino groups of lysine, arginine, and histidine side chains (T.E. Creighton, *Proteins: Structure and Molecular Properties*, W.H. Freeman & Co., San Francisco, pp. 79-86 [1983]), acetylation of the N-terminal amine, and amidation of any C-terminal carboxyl group.
  - [0036] The availability of a DNA sequence encoding M-CSF permits the use of various expression systems to produce the desired polypeptides. Construction of expression vectors and recombinant production from the appropriate DNA sequences are performed by methods well known in the art. These techniques and various other techniques are generally performed according to Sambrook et al., Molecular Cloning A Laboratory Manual, Cold Spring Harbor, New York (1989), and Kriegler, M., Gene Transfer and Expression, A Laboratory Manual, Stockton Press, New York (1990).

#### Purification of M-CSF

[0037] Purification steps are employed to ensure that the M-CSF is isolated, prior to crystallization, in a relatively homogeneous state. In general, a higher purity solution increases the likelihood of success of subsequent crystallization steps. Typical purification methods include the use of centrifugation, partial fractionation using salt or organic compounds, dialysis, conventional column chromatography (such as ion-exchange, molecular sizing chromatography etc.), high performance liquid chromatography (HPLC), and gel electrophoresis methods (see, e.g., Deutcher, "Guide to Protein Purification" in Methods in Enzymology (1990), Academic Press, Berkely, CA, which is incorporated herein by reference for all purposes). Preferred purification conditions for generating unusually homogeneous M-CSF species as well as purification of these species are disclosed, for example, in U.S. Patent No. 4,929 700. Other purification methods are known and will be apparent to one of skill in the art.

## 55 Crystallization of M-CSF

[0038] Although many of the same physical principles govern crystallization of polypeptides (including M-CSF dimers) and small molecules, the actual crystallization mechanisms differ significantly. For example, the lattice of small-

molecule crystals effectively excludes solvent while that of polypeptide crystals includes substantial numbers of solvent molecules. Thus, special techniques must typically be applied to crystallize polypeptides.

[0039] Polypeptide crystallization occurs in solutions where the polypeptide concentration exceeds its solubility maximum (i.e., the polypeptide solution is supersaturated). Such "thermodynamically metastable" solutions may be restored to equilibrium by reducing the polypeptide concentration, preferably through precipitation of the polypeptide crystals. Often polypeptides may be induced to crystallize from supersaturated solutions by adding agents that alter the polypeptide surface charges or perturb the interactions between the polypeptide and bulk water to promote associations that lead to crystallization.

[0040] Compounds known as "precipitants" are often used to decrease the solubility of the polypeptide in a concentrated solution. Precipitants induce crystallization by forming an energetically unfavorable precipitant-depleted layer around the polypeptide molecules. To minimize the relative amount of this depletion layer, the polypeptides form associations and ultimately crystals as explained in Weber, *Advances in Protein Chemistry* 41:1-36 (1991). In addition to precipitants, other materials are sometimes added to the polypeptide crystallization solution. These include buffers to adjust the pH of the solution (and hence-surface charge on the peptide) and salts to reduce the solubility of the polypeptide. Various precipitants are known in the art and include the following: ethanol, 3-ethyl-2,4 pentanediol; and many of the polyglycols, such as polyethylene glycol. A suitable precipitant for crystallizing M-CSF is polyethylene glycol (PEG), which combines some of the characteristics of the salts and other organic precipitants (see, for example, Ward et al., J. Mol. Biol. 98:161 [1975] and McPherson J. Biol. Chem. 251:6300 [1976],).

[0041] Commonly used polypeptide crystallization methods include the following techniques: batch, hanging drop, seed initiation, and dialysis. In each of these methods, it is important to promote continued crystallization after nucleation by maintaining a supersaturated solution. In the batch method, polypeptide is mixed with precipitants to achieve supersaturation, the vessel is sealed and set aside until crystals appear. In the dialysis method, polypeptide is retained in a sealed dialysis membrane which is placed into a solution containing precipitant. Equilibration across the membrane increases the polypeptide and precipitant concentrations thereby causing the polypeptide to reach supersaturation levels.

20

25

[0042] In the hanging drop technique, an initial polypeptide mixture is created by adding a precipitant to concentrated polypeptide solution. The concentrations of the polypeptide and precipitants are such that in this initial form, the polypeptide does not crystallize. A small drop of this mixture is placed on a glass slide which is inverted and suspended over a reservoir of a second solution. The system is then sealed. Typically the second solution contains a higher concentration of precipitant or other dehydrating agent. The difference in the precipitant concentrations causes the protein solution to have a higher vapor pressure than the solution. Since the system containing the two solutions is sealed, an equilibrium is established, and water from the polypeptide mixture transfers to the second solution. This equilibration increases the polypeptide and precipitant concentration in the polypeptide solution. At the critical concentration of polypeptide and precipitant, a crystal of the polypeptide will form. The hanging drop method is well known in the art (see, McPherson J. Biol. Chem. 251:6300 [1976],).

[0043] Another method of crystallization introduces a nucleation site into a concentrated polypeptide solution. Generally, a concentrated polypeptide solution is prepared and a seed crystal of the polypeptide is introduced into this solution. If the concentrations of the polypeptide and any precipitants are correct, the seed crystal will provide a nucleation site around which a larger crystal forms.

[0044] In preferred embodiments, the crystals of the present invention will be formed from a dimer of M-CSF polypeptides. Preferred crystals are typically at least about 0.2 x 0.2 x 0.05 mm, more preferably larger than 0.4 x 0.4 x 0.4 mm, and most preferably larger than 0.5 x 0.5 x 0.5 mm. After crystallization, the protein may be separated from the crystallization mixture by standard techniques.

[0045] The crystals so produced have a wide range of uses. For example, high quality crystals are suitable for X-ray or neutron diffraction analysis to determine the three-dimensional structure of the M-CSF and, in particular, to assist in the identification of its receptor binding site. Knowledge of the binding site region and solvent-accessible residues available for contact with the M-CSF receptor allows rational design and construction of agonists and antagonist for M-CSFs. Crystallization and structural determination of M-CSF muteins having altered receptor binding ability or bioactivity allows the evaluation of whether such changes are caused by general structural deformation or by side chain alteration at the substitution site.

[0046] In addition, crystallization itself can be used as purification method. In some instances, a polypeptide or protein will crystallize from a heterogeneous mixture into crystals. Isolation of such crystals by filtration, centrifugation, etc. followed by redissolving the polypeptide affords a purified solution suitable for use in growing the high-quality crystals necessary for diffraction studies. These high-quality crystals may also be dissolved in water and then formulated to provide an aqueous M-CSF solution having various uses known in the art including pharmaceutical purposes.

[0047] Of course, amino acid sequence variants of M-CSF may also be crystallized and used. These mutants can be used for, among other purposes, obtaining structural information useful for directing modification of the binding affinity for M-CSF receptors. As with the naturally occurring forms, the modified M-CSF forms may be useful as pharmaceu-

tical agents for stimulating bone marrow proliferation, overcoming immune suppression and fungal diseases induced by chemotherapy, improving therapeutic efficacy, and lessening the severity or occurrence of side effects during therapeutic use of the present invention. Furthermore, modified M-CSFs may be useful for treatment of disease in which soluble or membrane-bound M-CSF causes or exacerbates the disease state.

#### Characterization of M-CSF

5

[0048] After purification, crystallization and isolation, the subject crystals may be analyzed by techniques known in the art. Typical analysis yield structural, physical, and mechanistic information about the peptides. As discussed above, X-ray crystallography provides detailed structural information which may be used in conjunction with widely available molecular modeling programs to arrive at the three-dimensional arrangement of atoms in the peptide. Exemplary modeling programs include "Homology" by Biosym (San Diego, California), "Biograf" by BioDesign, "Nemesis" by Oxford Molecular, "SYBYL" and "Composer" by Tripos Associates, "CHARM" by Polygen (Waltham, Massachusetts), "AMBER" by University of California, San Francisco, and "MM2" and "MMP2" by Molecular Design, Ltd.

[0049] Peptide modeling can be used to design a variety of agents capable of modifying the activity of the subject peptide. For example, using the three-dimensional structure of the active site, agonists and antagonists having complementary structures can be designed to enhance the therapeutic utility of M-CSF treatment or to block the biological activity of M-CSF. Further, M-CSF structural information is useful for directing design of proteinaceous or non-proteinaceous M-CSF agonists and antagonists, based on knowledge of the contact residues between the M-CSF ligand and its receptor. Such residues are identified by the M-CSF crystal structure as those which are solvent-accessible, distal to the carboxyl terminal membrane anchoring region not involved in dimer interface stabilizations, and possibly including residues not conserved between human and mouse M-CSF (which does not recognize the human M-CSF receptor).

### Example 1

[0050] Systematic crystallization trials with M-CSF were made using the hanging drop technique. A microdroplet (5µl) of mother liquor is suspended from the underside of a microscope cover slip, which is placed over a well containing 1ml of the precipitating solution. 60-70 initial trials were set up, in which pH, temperature, counterion and precipitant were varied. From these trials, the few that gave promising microcrystals were picked for more careful examination.

[0051] It was discovered that suitable crystals may be grown from a 20  $\mu$ l drop containing: 10 mg/ml protein, 100 mM MgCl<sub>2</sub>, 50 mM Tris.Cl, pH 8.5, and 12% PEG 4000. This drop was equilibrated against a reservoir containing 24% PEG 4000. Tiny, needle-like crystals appeared in 2-3 days which were redissolved in 10  $\mu$ l water and recrystallized at room temperature. Good quality chunky crystals appeared in 7-9 days in sizes ranging from 0.3 x 0.3 x 0.3 mm to 0.5 x 0.5 x 1.0 mm.

[0052] Precession photographs revealed the space group to be  $P2_12_12_1$  with unit cell dimensions: a=33.53Å, b=65.11Å, c=159.77Å. This gives a unit cell volume of 349084.5Å<sup>3</sup>, which is consistent with a dimer in the crystallographic asymmetric unit, and 52% of the unit cell volume being occupied by solvent. The crystals diffracted to a resolution of 3Å on a Rigaku rotating anode X-ray generator (Danvers, Massachusetts) operated at 50 kV and 60mA, and to 2.6Å in synchrotron radiation.

[0053] Screening for heavy atom derivatives was done by soaking crystals into solutions of heavy-metal salts. Zero-level precession pictures of the soaks were used to identify potential derivatives. About 30 different soaking conditions were examined, of which 4 potential derivatives were identified. Unfortunately, some soaks caused the crystals to exhibit non-isomorphism (i.e., the heavy atom soaks induced a change in cell dimensions, making them unusable for phase calculation).

[0054] Three-dimensional intensity data were collected on film using an oscillation camera on the X-ray beam-line at the National Synchrotron Light Source, Brookhaven. Several other data sets, of native (underivatized M-CSF) as well as potential derivative crystals have been collected on a Rigaku X-ray generator. The following data sets were collected.

35

Crystal	Resolution(Å)	N(observations)	N(unique)	X-Ray Source
Native	2.8	27922	7311	Synchrontron
Native	2.9	35236	7002	Rigaku (film)
Native	3.5	5144	5116	Rigaku (diffractometer)

#### (continued)

Crystal	Resolution(Å)	N(observations)	N(unique)	X-Ray Source
K <sub>2</sub> Hg(SCN) <sub>4</sub>	3.5	15885	4119	Rigaku (film)
UO <sub>2</sub> Cl <sub>2</sub>	3.5	25492	5048	Rigaku (film)
Cis-Pd(NH <sub>3</sub> ) <sub>2</sub> Cl <sub>2</sub>	3.1	26122	6304	Synchrotron

#### Example 2

[0055]

5

10

[0055] Recombinant M-CSF polypeptides were purified from *E. coli.* and renatured to form a disulfide-linked dimeric protein as described in U.S. Patent No. 4,929,700. Crystallization of the resulting unglycosylated M-CSFα protein (amino acids 4-158 in homodimeric form) was performed by the hanging drop method. Glass microscope plates were siliconized prior to use by dipping immersion into a 1% (volume:volume) solution of the organosilane compound, Prosil-28 (PCR Incorporated, Gainesville, Florida, 32602) washing the treated glass plates with water, and baking at 180 degrees.

[0056] A 2 mg/ml aqueous solution of purified human recombinant M-CSF was dialyzed and concentrated against 50 mM Tris-HCl (pH 8.5) using a dialysis tubing having a 10kD cutoff. The final concentration of polypeptide (10 mg/ml) was determined by ultraviolet spectrophotometry at 280 nm.

[0057] About 7 microliters of the concentrated solution was mixed in each well of the spot plate with 7 microliter of 20% (v/v) PEG 4000, 0.2 M MgCl<sub>2</sub>, 0.1 M Tris-HCl (pH 8.5). The spot plate was then placed in a clear plastic sandwich box containing 20 ml of 23% PEG 4000, 0.2 M MgCl<sub>2</sub>, 0.1 M Tris-HCl (pH 8.5) and the box was immediately sealed and stored at room temperature. Minor variations in this procedure such as altering buffer conditions are within the scope of the present invention. For example, in a preferred embodiment of the present invention, buffer conditions were altered to include 150 mM MgCl<sub>2</sub> and 24% PEG 4000.

[0058] After 3-5 days, small microcrystals having a size of 0.1 x 0.1 x 0.05 mm appeared in each well. These microcrystals were isolated and redissolved in 25 microliter of 50 mM Tris-HCl and allowed to stand at room temperature. The purified M-CSF crystallized from solution into large hexagonal prism shaped crystals ranging in size from 0.3 x 0.3 x 0.3 mm to 1 mm x 2 mm x 0.5 mm. These crystals were stable at room temperature for at least three months. In some instances, an artificial mother liquor was prepared using 23% PEG 4000 and 150 mM MgCl<sub>2</sub> crystals were then added to this mother liquor. In these cases, the crystals were removed from the mother liquor immediately prior to analysis.

[0059] Using reducing and non-reducing SDS-PAGE analysis, the M-CSF in the crystals was shown to be identical in molecular weight to the biologically active starting material. Thus, the M-CSF structure obtained from the crystals is likely to be essentially identical to the structure of biologically active M-CSF.

#### Example 3

35

[0060] Glass microscope slides were prepared as described in Example 2. 7 microliters of the same concentrated M-CSFα protein solution was mixed in each well of the spot plate with 7 microliter of 30% (v/v) PEG 4000, 0.2 M ammonium acetate, 0.1 M acetate buffer (pH 7.5). The spot plate was then placed in a clear plastic sandwich box containing 20 ml of 30% PEG 4000, 0.2 M ammonium acetate, 0.1 M acetate buffer (pH 7.5) and the box was immediately sealed and stored at room temperature. After 3-5 days, thin, plate-like, fragile crystals having a size of approximately 0.3 x 0.3 x 0.05 mm appeared.

#### 45 Example 4

## Preliminary X-ray Analysis

[0061] X-ray crystallographic analysis using precession photographs showed that the crystals produced in Example 2 have an orthorhombic crystal lattice in the P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub> space group with cell dimensions a=33.54, b=65.26, c=159.63 d = 90.0, c = 90.0 and f = 90.0 angstroms and diffract to a nominal resolution of 2.6 angstroms using synchrotron radiation. These data provided a unit cell volume of 348084.5 angstroms<sup>3</sup>, which is consistent with a dimer in the crystallographic asymmetric unit with 52% of the unit cell volume being occupied by solvent. Fig. 1 is a 12-degree precession photograph of the Okl-zone section of the M-CSF crystal. The photograph was taken using a precession camera manufactured by Enraf-Nonius Company (Delft, Holland), mounted on a Rigaku RU-200 X-ray generator operated at 50kV and 50mA.

#### Example 5

# Testing of M-CSF Receptor Binding Ability Using Soluble M-CSF Receptor

An essential step in the biological function of M-CSF in vivo is the binding to the M-CSF receptor, also [0062]referred to as the c-ims gene product. Recombinant human soluble M-CSF receptor (rhsM-CSFR), representing amino acids 20 to 511 (Coussens, Let al, Nature, 320:277 (1986)) was used as an in vitro assay reagent to test the receptorbinding ability of M-CSF proteins. To generate a soluble form of the transmembrane receptor, only the extracellular domain of the human M-CSF receptor was expressed in a baculovirus/insect cell recombinant expression system. In order to purify the soluble receptor without adversely affecting tertiary or quaternary structure, non-denaturing chromatographic methods were chosen, as described below. Other choices exist for the purification of the recombinant receptor. Affinity chromatography may be employed when either a suitable antibody to or ligand for the receptor are available. Alternatively, "tags" may be added to the C-terminus of the recombinant receptor, i.e., KT3 antibody recognition sequence, and purified by an anti-tag antibody, i.e., KT3, column, for use in affinity chromatography. In expression systems in which the rhsM-CSFR is glycosylated, lectin chromatography can be used to enrich for specific glycoproteins. The rhsM-CSFR can be used to study ligand/receptor interactions as well as ligand-induced receptor dimerization. The assay used to detect ligand/receptor binding employed the use of size exclusion-HPLC, essentially as described in European Patent Application WO92/21029, C. Cunningham, et al., with the following modifications: the column used was a Superose 6 (Pharmacia LKB Biotechnology, Inc.) and the mobile phase was PBS at 0.5 ml/min and a M-CSF to rhsM-CSFR ratio of 1:2. At this ratio, the M-CSF/rhsM-CSFR complex chromatographed with an apparent hydrodynamic radius of 190,000 molecular weight, the molecular weight expected for a M-CSF(rhsM-CSFR)2 complex. Other assays may be employed to measure ligand/receptor binding or receptor dimerization such as chemical crosslinking and SDS-PAGE or immunoprecipitation and SDS-PAGE. Molecules that inhibit receptor dimerization but not ligand binding provide another method to antagonize M-CSF actions.

[0064] The DNA encoding rhsM-CSFR was cloned for expression in insect cells using the following general strategy. The portion of the *c-fms* gene corresponding to amino acids one to 511 was amplified from human macrophage cDNA by polymerase chain reaction (PCR) using an upstream primer of: 5'-GCGTACCATGGGCCCAGGAGTTCTGC-3' (SEQ ID NO.9) and a downstream primer of: 5'-AGTCGAGGATCCTCAATCCGGGGGATGCGTGTG-3' (SEQ ID NO. 10). The underlined sequences are the Ncol and BamHI restriction sites used to subclone the PCR product into the pAcC5 vector (Luckov et al, Bio/Technology 6:47-55). The pAcC5:hsM-CSFR vector was expressed in Sf9 insect cells using a baculovirus helper vector as previously described (Summers, et al., A Manual of Methods for Baculovirus Vectors and Insect Cell Culture Procedures (1987)).

Approximately two liters of serum-free 72-hour conditioned medium was collected by centrifugation and fil-[0065]tration from Sf9 cells infected with pAcC5:hsM-CSFR construct described above. The material was diafiltered with DEAE buffer A [10 mM Tris, pH 8.8, containing the following protease inhibitors (which were added all buffers throughout the purification): 1 mM EDTA, 2 μg/ml leupeptin and 100 μM PMSF] and concentrated 20-fold with a 20,000 molecular weight cut-off Pyrostat Ultrafiltration Membrane (Sartorius). The retentate was loaded onto a DEAE Sepharose column (Pharmacia LKB Biotechnology, Inc., Piscataway, New Jersey) having a bed volume of 100 ml that had been pre-equilibrated with DEAE Buffer A. Elution was at 5 ml/min with a 0-0.8 M NaCl gradient in 500 ml of DEAE Buffer A. Fractions enriched in rhsM-CSFR were detected by Western Analysis [Burnett, R., Anal. Biochem., 112:195 (1981)] and dot blot analysis of serially diluted fractions, using anti-c-FMS monoclonal antibodies (Oncogene Sciences, Inc.). The dot blot assay was used throughout the purification to identify fractions containing rhsM-CSF. Enriched fractions were pooled, made 0.8 M in ammonium sulfate, adjusted to pH 7.0 and loaded onto a Phenyl TSK-5-PW HPLC column (7.5 x 75 mm) (BioRad). The column was eluted at 1 ml/min with a decreasing ammonium sulfate gradient over 45 minutes, peak fractions were pooled and concentrated 10-fold with a stir cell concentrator using a YM30 membrane (Amicon). The retentate was chromatographed with FG3000XL size exclusion column (DU PONT, Wilmington, Delaware) using a mobile phase a phosphate-buffered saline (PBS) at 3 ml/min. The purified receptor was pooled, concentrated to 1 mg/ml as above and stored at 4°C. This process recovered 650 µg of rhsM-CSFR, purified 200-fold. The preparation was about 95% homogeneous as assayed by SDS-PAGE stained with Coomassie Blue.

#### Example 6

50

# Crystallization of M-CSF/Soluble M-CSF Receptor Complex

To crystallize the M-CSF/rhsM-CSFR complexes, glass microscope slides are prepared as described in Example 2. The M-CSF composition used is incubated with a purified soluble form of the M-CSF receptor, truncated at a residue before the transmembrane region, to form an M-CSF/ receptor complex. In certain cases, the rhsM-CSFR is deglycosylated prior to the size exclusion step by incubation with N-glycanase (Genzyme, Cambridge Massachusetts)

according to the manufacturer's instructions. A small quantity of M-CSF/receptor solution is mixed in each well of the spot plate with a comparable quantity of a drop solution (such as about 20% (v/v) PEG 4000, 0.2 M MgCl<sub>2</sub>, 0.1 M Tris-HCl (pH 8.5)). The spot plate is then placed in a clear plastic sandwich box containing a small amount of precipitant solution (such as about 23% PEG 4000, 0.2 M MgCl<sub>2</sub>, 0.1 M Tris-HCl (pH 8.5)). The box is immediately sealed and stored at room temperature.

[0067] After a few days, crystalline M-CSF-receptor complex is isolated and redissolved in a solution containing about 50 mM Tris-HCl and is allowed to stand at room temperature. The purified M-CSF-receptor complex crystallizes from solution to form crystals for X-ray structural analysis. To facilitate solution of the crystal structure of such complexes, truncated, non-glycosylated forms of the rhsM-CSFR (described above) which retain M-CSF binding ability may be employed to generate M-CSF-receptor complex crystals.

#### Example 7

[0068] The biological activity of the non-glycosylated, truncated sequence used in Examples 2 and 3 was shown to be equal to that of the mature protein purified from human urine (Halenbeck, R., et al., Bio/Technology, 7:710-715 [1989]). As noted, the resulting crystals had an orthorhombic crystal lattice in the P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub> space group, with cell dimensions a=33.54, b=65.26, and c=159.63 Å. Intensity data were collected using imaging plates mounted on a Weissenberg camera modified for macromolecular crystallography at the Photon Factory in Tsukuba, Japan. Native data to a nominal resolution of 2.0Å, and mercury and platinum derivative data were collected using 1.0Å radiation. Two crystal settings were used to collect native data [Rmerge (I)=7.0%, using all measurements with I>0.0].

[0069] Heavy atom derivatives of M-CSF crystals were prepared by soaking crystals in heavy atom compounds dissolved in the reservoir solution. Isomorphous and anomalous difference Patterson maps clearly revealed one site for the mercury and two sites for the platinum derivative. Anomalous and isomorphous phase information as used in initial phase refinement with the PROTEIN program package. The final figure of merit was 0.62 (8.0-3.0Å, 6960 reflections). After solvent flattening, B.C. Wang, *Methods Enzymol*, 115:90 (1985), two bundles of four alpha helices related by an approximate two-fold axis could be seen in the electron density map. Rotation and translation parameters of this non-crystallographic axis were refined by a density correlation method, J.M. Cox, *J. Mol. Biol.* 28:151 (1967). Phases were then iteratively refined by molecular avenging and solvent flattening, G. Bricogne, *Acta Cryst.*, 32:832 (1976), using an envelope calculated by putting 5Å spheres around all the atoms in the four helical bundle. Chain tracing and model building were done in the resulting map, (using the program FRODO), T.A. Jones, *Methods Enzymol.* 115:157 (1985), keeping the original MIR map as a reference.

[0070] The starting partial model for refinement contained only a polyalanine backbone for eight helices making up the two bundles. Positional refinement using the program XPLOR, A.T. Brunger, *J. Mol. Biol.* 203:803 (1985), gave an R-factor of 0.49 to 3.0Å. Phase combination with the refined MIR phases resulted in a map of sufficient quality to allow the tracing of two long loops traversing the four helical bundle and a short loop connecting two of the helices. Two strong peaks in the density, one at the top of the first helix, and the second lying directly on the molecular two fold axis, were assigned as disulfide bonded cysteines. The number of residues between these two peaks uniquely identified the position in the sequence of these cysteines and consequently the sequence of the intervening residues. This initial registration was confirmed by the presence of a number of regions of strong density corresponding to aromatic side chains in the sequence. Partial model phase combination using the added loops and those side chains that were visible allowed the remaining residues to be registered, thus determining the overall topology of the molecule. The presence of seven disulfide bonds in the dimer served as important "tether points" to confirm the correctness of the tracing.

[0071] As shown in Fig. 2, the overall topology of this form of M-CSF is that of an antiparallel four  $\alpha$ -helical bundle, in which the helices run up-up-down-down, unlike the more commonly observed up-down-up-down connectivity of most four helical bundles. A long crossover connection links helix A to helix B and a similar connection is found between helices C and D.

[0072] A striking difference from other cytokines and other four helix bundle structures is that the truncated M-CSF $\alpha$  forms a disulfide-linked dimer, in which the bundles are linked end-to-end, forming an extremely flat, elongated structure (approximate dimensions 85 x 35 x 25 Å) as shown in Figs. 3 and 4A and 4B. There are three intramolecular disulfide bonds in each monomer (Cys7-Cys90, Cys48-Cys139, Cys102-Cys146) all of which are at the distal end of the molecule. One interchain disulfide bond (Cys31-Cys31) is located at the dimer interface with the noncrystallographic two-fold symmetry axis passing through it as shown in Figs. 3 and 4A and 4B. Mutation experiments indicate that all of the cysteine residues in this form of M-CSF may be necessary for full biological activity. The structure described herein suggests that their role is primarily structural rather than being related to receptor recognition.

[0073] Appendix 1 provides the three-dimensional structure of the truncated recombinant M-CSFα dimer as identified by the alpha-carbon positions of the amino acid residues in the sequence. The five carboxy terminal amino acids of each polypeptide of the dimer were not included. As will be recognized to those of skill in the art, the information in Appendix 1 is provided in the format used by the Brookhaven Protein Data Bank.

[0074] As shown, the molecule has an unusual topology which identifies important regions of M-CSF with regard to M-CSF receptor binding. Specific residues in helices A, C, and D appear to be involved in the specificity of the interaction. Altering solvent accessible residues in these regions by site directed mutagenesis to increase or decrease side-chain interactions with the receptor may be useful to generate M-CSF agonists or antagonists. For example, changing one or more histidines to non-hydrogen-donor amino acids of similar size may create an M-CSF with altered receptor binding ability.

#### Example 8

25

# Preparation of M-CSF Heterodimers

# Purification of M-CSF Monomers:

[0075] E. coli harboring the pL-M-CSF vector for NΔ3CΔ158M-CSFα described in U.S. Patent No. 4,929,700 or NΔ3CΔ221 M-CSFβ C157S, C159S (Kawasaki, et al., in Colony-Stimulating Factors, Dexter, T., Garland, J., and Testa, N., eds. [1990]) were grown in 1 liter of minimal salts medium containing glucose and the appropriate antibiotic. Expression of M-CSF was induced by shifting the temperature to 39°C for 4 hr. following the addition of casamino acids to 2%. The cells were harvested by centrifugation and lysed by sonication in 50 mM Tris (pH 8.5), 10 mM EDTA. The cell debris paste was solubilized in 8 M urea under reducing conditions. After incubation at 37°C for 30 min., the solubilized M-CSF was clarified by centrifugation and filtration, and then loaded onto a Bio-Gel TSK-DEAE-5-PW column (7.5 x 75 mm) (BioRad Laboratories, Richmond, California) equilibrated in 8 M urea in 10 mM Tris (pH 8.5), 5 mM DTT, 0.5 mM EDTA. The monomeric M-CSF was eluted with a 45-min, 0-0.6 M NaCl gradient. The M-CSF peak fractions were pooled and concentrated to 10 mg/ml with a Centricon 10 microconcentrator (Amicon).

# Formation and Analysis of Active M-CSF Heterodimers

[0076] The M-CSF homodimers were refolded by diluting to a protein concentration of 0.5 mg/ml in precooled 50 mM Tris (pH 8.5), 5 mM EDTA, 2 mM reduced glutathione and 1 mM oxidized glutathione, and then incubating at 4°C. The heterodimer was refolded by diluting 158 and 221F monomer pools to 1 mg/ml in the same buffer. To monitor the refolding, size exclusion high pressure liquid chromatography (SE-HPLC) analysis was performed by immediately injecting reaction samples onto a G3000SW Ultropack (Pharmacia LKB Biotechnology, Inc., Piscataway, New Jersey) column size-exclusion (7.5 x 600 mm) equilibrated in PBS (pH 6.8).

[0077] Fractionated products were analyzed on reducing SDS-PAGE and stained with Coomassie, according to the method of Laemmli, *Nature* (Canada) 227: 680-685 (1970). Biological activity was determined using the M-CSF dependent NFS-60 bioassay (see Example 10 below). Antibody neutralization experiments were carried out by preincubating approximately 5,000 units of M-CSF dimer with varying dilutions of the neutralizing M-CSF 5H410 Mab (made to refolded *E. coli* CΔ 150 M-CSFα dimer) prior to bioassay. (Halenbeck *et al.*, Bio/Technology 7:710 (1989)

[0078] The heterodimeric M-CSF product was designed to consist of one chain of short clone (from amino acid 4 to 158) and one chain of long clone (from amino acid 4 to 221). The long-clone chain (221F) also contained substitutions of serine for the two non-essential cysteines (at 157 and 159) to minimize the possibility of higher-order oligomer

[0079] Solubilized refractile bodies of M-CSF 158 and 221F were separately chromatographed by DEAE-HPLC in 8 M urea. Only one major protein peak eluted in each case, and the peak fractions were pooled, based on an analysis of purity by non-reducing SDS-PAGE and Coomassie staining (data not shown). The resulting monomer was over 90% pure in each case. The monomers were separately concentrated to 10 mg/ml, diluted in refolding buffer, and refolded at 4°C.

[0080] To compare the rates of dimerization of short- and long-clone M-CSF, 20 µl of each refolding reaction was injected on a SE-HPLC column at 0, 2, 18 and 72 hr. The amount of dimeric M-CSF formed was determined from the peak area at the molecular weight expected for dimer. In both refolding reactions the M-CSF was mostly equilibrated to monomer at t=0 and had become about 40% dimeric by 2 hr and nearly 75% dimeric by 18 hr. The similarity of the ratio of dimer to monomer between the refolded 158 and 221F strongly suggests that the rate of dimer formation is the same for long- and short-clone M-CSF. Thus, when equal moles of 158 and 221F are present in a refolding reaction, the final relative ratios of 158 homodimer to 221F homodimer to 158/221F heterodimer are predicted to be 1:1:2 (Similar distributions have been observed *in vivo* for isozymes of lactate dehydrogenase.)

## Biological Activity of Refolded Homodimers and Heterodimers

[0081] The biological activity of the refolded homodimers and heterodimers described above, was examined using the *in vitro* M-CSF-dependent NFS-60 bioassay (see Example 10 below). Figures 5A-5C show the result of these studies. These SE-HPLC and biological activity profiles analyzed after 72 hr of refolding, show that the heterodimer, Figure 5C, displays activity very similar to that of the two homodimers, (Figures 5A and 5B). Given that the separation of the heterodimer from the homodimers was nearly complete, it can be concluded that the heterodimer is fully biologically active *in vitro*.

[0082] To verify that the M-CSF protein eluting from these columns at the predicted heterodimeric position (between the two homodimers) actually did consist of equal moles of short- and long-clone monomers, analysis of a preparative purification of the 158/221F heterodimer was carried out. Phenyl-HPLC was performed as described above and was shown to completely resolve the heterodimer from the 158 and 221F homodimers, as seen in Figure 6A.

#### Preparative Purification of M-CSF Heterodimers:

[0083] The refolded M-CSF was adjusted to pH 7.0 with 1 N HCl, and ammonium sulfate was added to 1.2 M. The protein was loaded onto a Bio-Gel TSK-Phenyl-5-PW column (7.5 x 75 mm) (BioRad, Richmond, California) equilibrated in 1.2 M ammonium sulfate, 100 mM phosphate (pH 7.0). The M-CSF was eluted with a decreasing gradient of ammonium sulfate from 40% to 80% buffer B (10 mM phosphate, pH 7.0) in 24 min.

[0084] Reducing and non-reducing SDS-PAGE (Figures 6B and C) showed that internal controls (the 158 and 221F dimers) were purified to approximately 95% homogeneity by this column, and each consisted of the single expected monomeric band. The gel analysis also showed that the heterodimer was purified to approximately 95% homogeneity and that it consisted of equivalent amounts of 158 and 221F monomers. Recovery of purified 158/221F heterodimer from refractile body paste to final product was greater than 15%.

**[0085]** The bioactivity of the dimeric M-CSF species was determined and, when compared to the  $A_{280}$  profile in Figure 6A, confirms the finding that the heterodimer is fully active. The specific activity of the 158/221F heterodimer, calculated using the peak fraction, was  $8.0 \times 10^7$  units/mg, compared to  $9.0 \times 10^7$  and  $6.8 \times 10^7$  units/mg for 158 and 221F homodimers, respectively.

[0086] The biological activity of all three dimer species was neutralized to the same extent in serial dilution neutralization experiments using the 5H410 M-CSF Mab in the NFS-60 bioassay. This antibody also neutralizes "naturally refolded" Chinese hamster ovary cell (CHO)-expressed M-CSF in a similar fashion. This result further suggests that the refolded conformation of the new M-CSF heterodimer is essentially native-like, at least with regard to the region within the first 150 amino acids that is responsible for *in vitro* activity.

#### 35 Example 9

15

25

50

# Selection of Amino Acid Substitutions in M-CSF Based on Crystallographic Data

[0087] The X-ray crystallographic data described above provided sufficient structural information regarding M-CSF to be able to identify a limited subset of the amino acids in the protein that are likely to be crucial for M-CSF receptor binding and biological activity and thus which represented likely candidates for mutagenesis with the ultimate goal of providing M-CSF muteins having altered biological activity (*i.e.*, agonists or antagonists). Based on this information, several criteria were used to generate a list of possible target amino acids for substitution.

[0088] The first criterion was solvent exposure or solvent accessibility, which refers to amino acids residues at the surface of the protein. Residues having a solvent accessible surface area of greater than about 0.25 and preferably greater than about 0.4 are preferred based on normalization of the surface area of the amino acid accessible when in the trypeptide gly-x-gly (Kabsch, W. et al., Biopolymers 22:2577 (1983)). Residues were chosen which do not interact with other parts of the protein such as the dimer interface in order to maintain the relative orientation of monomers and to avoid disturbing the process of protein folding. Still another criterion used in certain instances in selecting candidate amino acid substitutions is the relationship of the residues to corresponding residues in mouse M-CSF. Another important selection criterion was that the substitutions be non-conservative so as to attempt to disrupt possible hydrogen bonding or hydrophobic interactions with M-CSF receptor residues.

[0089] Table 1 lists exemplary amino acid residues and exemplary substitutions. Using the criteria for selecting candidates for substitutions set forth above, those of ordinary skill in the art may readily ascertain other possible candidates for substitution.

Table I

Candidate Substitutions							
	Amino Acid and cation	Substitutions	٦				
	cation						
His (H)	15	Ala(A) or Leu(L)	7				
Gln (Q)	17	Ala(A) or Glu(E)					
Gln (Q)	79	Ala(A) or Asp(D)					
Arg (R)	.86	Glu(E) or Asp(D)					
Glu (E)	115	Ala(A)					
Glu (E)	41	Lys(K) or Arg(R)					
Lys (K)	93	Ala(A) or Glu(E)					
Asp (D)	99	Lys(K) or Arg(R)					
Leu (L)	55	Gln(Q) or Asp(N)					
Ser (S)	18	Ala(A) or Lys(K)	l				
Gin (Q)	20	Ala(A) or Asp(D)					
Arg (R)	21	Ala(A), Glu(E), or Asp(D)					
lle (I)	75	Lys(K) or Glu(E)					
Val (V)	78	Lys(K) or Arg(R)					
Leu (L)	85	Glu(E) or Asn(N)					
Asp (D)	69	Lys(K) or Arg(R)					
Asn (N)	70	Ala(A) or Glu(E)					
His (H)	9	Ala(A) or Asp(D)					
Asn (N)	63	Lys(K) or Arg(R)					
Thr (T)	34	Gin(Q) or Lys(K)					

[0090] It is not expected that every candidate substitution listed will result in the production of M-CSF agonists or antagonists (see Example 12 below). Rather they represent a non-exclusive list of candidates likely to result in the production of agonists or antagonists based on the selection criteria set forth above. It should also be noted that even if a variant does not act as an agonist or antagonist when compared with native M-CSF, the variant is still useful for conventional uses of the ligand (if it retains the same activity as the ligand) or as for example, a diagnostic reagent.

#### Example 10

5

10

15

20

25

30

35

45

# Preparation of H9A, H15A M-CSF Muteins

[0091] A variety of M-CSF muteins with altered solvent-accessible residues from regions of the M-CSF mature N terminus and helices A, C, and D were constructed using techniques known in the an. For example, two histidines in the N-terminal/A helix region were changed to alanine through site-directed mutagenesis of a truncated form of M-CSFα (encoded by pLCSF158A). Involvement of one of three M-CSF histidine residues in M-CSF receptor interaction was implicated by our observation that diethylpyrocarbonate (DEPC) modification of histidines in M-CSF at a 1:100 DEPC:histidine ratio (as described in *Meth. in Enzymol.* 47:431 (1977)) significantly reduced bioactivity.

[0092] Plasmid DNA pLCSF158A was prepared from the *E. coli* strain HW22 carrying the plasmid pLCSF158A

[U.S. Patent No. 4,929,700, Example 6, "E. coli strain HW22 transformed with pJN653 containing the asp<sub>59</sub>SCSF/NΔ3CΔ158 gene"). The strain was grown in 350 ml R2 media (2X Luria Broth containing 1% sodium chloride and no glucose, J. Bact., 74:461 (1957)) containing 50 micrograms/ml ampicillin at 30° C with shaking overnight. Plasmid DNA was prepared from the cells using a Qiagen-tip 100 column according to the manufacturer's directions.

[0093] Twenty micrograms of pLCSF158A DNA were digested with 66 units of *Hind*III and 66 units of *Stu*I at 37°C for 3hr. 20 min. in 200 microliters 1X New England Biolabs NEBuffer 2 (New England Biolabs, Beverly, Massachusetts). The DNA was extracted with phenol and chloroform, then ethanol precipitated. The DNA was treated with one unit of Calf Intestinal Alkaline Phosphatase in 100 microliters of 1X Dephosphorylation Buffer, supplied by Boehringer Mannheim (Indianapolis, Indiana), at 37°C for 30 min. An additional unit of Calf Intestinal Alkaline Phosphatase was added to the reaction and incubation was continued at 50°C for 1 hr.

[0094] The resulting DNA was then run on a 1 % FMC Bioproducts (Rockland, Maine) Sea KEM<sup>®</sup> GTC<sup>®</sup> agarose gel. The 5.7 kb pLCSF158A fragment was cut from the gel and purified on Qiagen (Chatsworth, California) Qiaex beads according to the manufacturer's directions.

[0095] Polymerase chain reaction (PCR) was then performed and a PCR product was produced that contained a mutagenized M-CSF sequence in which histidines 9 and 15 (counting from the mature N-terminus) were altered to alanine (generating an H9A, H15A PCR fragment). The 5' portion of the M-CSF gene was amplified from the plasmid pLCSF158A in a PCR reaction using the primers LF73 and LF74. Details of PCR are provided by Mullis, K. et al., U.S. Patent No. 4,683,202; Ehrlich, H., U.S. Patent No. 4,582,788; Saiki et al., U.S. Patent No. 4,683,195; Mullis, K.B., Cold Spring Harbor Symp. Quant. Biol. 51:263-273 (1986); Saiki et al., Bio/Technology 3:1008-10012 (1985); and Mullis, K.B. et al., Meth. Enzymol 155:335-350(1987), all of which are incorporated herein by reference. The sequences of these primers are:

## 20 LF73:

5'-AGGTGTCTCATAGAAAGTTCGGACGCAGGCCTTGTCATGCTCTTCATAATCCTTGG-3' (SEQ ID NO. 1)

25 LF74:

5'-CAGGAGAAAGCTTATGTCTGAATATTGTAGCGCCATGATTGGGAGTGGAGCCCTGCAG-3' (SEQ ID NO. 2)

30

[0096] The expected PCR product was designed to include 337 bp of pLCSF158A sequence, with the *Hind*III and *Stu*I sites located at each end of the product for cloning, and the histidine codons for His 9 and His 15, CAC, mutated to an alanine codons, GCC.

This product was amplified in four separate PCR reactions each containing 100 ng of pLCSF158A DNA, 50 pmoles LF73; 50 pmoles LF74, 37.5 µM dNTPS, 5% glycerol, 1X Perkin Elmer Cetus PCR Buffer, and 2.5 units of Perkin Elmer Cetus AmpliTaq<sup>®</sup> DNA Polymerase in a 100 microliter volume. The amplification was carried out in a Perkin Elmer Cetus DNA thermocycler. Before adding the AmpliTaq<sup>®</sup>, the reactions were brought to 95°C. The amplification was carried out for 25 cycles ramping to a denaturation temperature of 95°C in 1 sec., denaturing at 95°C for 1 min; ramping to an annealing temperature of 68°C in 1 sec., annealing at 68°C 1 min.; ramping to an extension temperature of 72°C in 30 sec., extending at 72°C 1 min., 30 sec. Final extension was carried out at 72°C for 10 min.

[0098] Five microliters of each reaction were run on a 3% agarose gel (1.5% FMC Bioproducts SeaKem<sup>®</sup> GTG<sup>®</sup> agarose, 1.5% FMC Bioproducts NuSeive<sup>®</sup> GTG<sup>®</sup> agarose in Tris-Borate buffer) (FMC Bioproducts, Rockland, Maine). Gels were then stained with ethidium bromide. For each reaction, a major band was visible at approximately 337 bp.

[0099] The four reactions were pooled, extracted with phenol and chloroform, precipitated with ethanol, resuspended and digested with 250 units of *Stul* in a final volume of 500 microliters 1X NEBuffer 2 at 37°C for 2hr., 500 units of *Hind*III were added to the reaction, the volume increased to 1 ml in 1X NEBuffer 2 and digestion was continued at 37°C for an additional 2.5hr. The DNA was electrophoresed on a 3% agarose gel. The 300 bp digested product was cut from the gel and purified on Qiagen Qiaex beads according to the manufacturer's directions.

[0100] Approximately 68 ng of the *HindIII/StuI* digested PCR product was then ligated to approximately 28 ng of the 5.7 kb *HindIII/StuI* digested pLCSF158A vector DNA at an insert-to-vector ratio of approximately 5:1. Ligation was carried out with 1 unit of Boehringer Mannheim T4 DNA ligase in 1X ligation buffer, supplied by the manufacturer, in a 20-microliter volume at 16°C overnight. As a control 28 ng of the 5.7 kb *HindIII/StuI* digested pLCSF158A vector DNA was ligated to itself under the same conditions with no insert present.

[0101] Half of each ligation mixture was used to transform competent E. coli DG116 (ATCC# 53606) cells using a protocol similar to the calcium chloride procedure described in Molecular Cloning a Laboratory Manual Maniatis et al., Cold Spring Harbor Laboratory (1982). Transformed cells were allowed to express at 30°C with no selection for 90 min., plated on R2-4 (10g tryptone, 5g yeast extract, 5g NaCl, 2 drops antifoam A, 4 ml 50% glucose and 15g agar in 1 liter)

plates containing 50 micrograms/ml ampicillin. The plates were incubated at room temperature 72 hr. One fourth of each transformation was plated. For the ligation containing the insert, 66 ampicillin resistant colonies appeared on the plates. For the ligation with no insert, no colonies appeared.

[0102] One of these colonies, designated strain TAF172-2, was picked and cultured in 350 ml R2 broth with 50 micrograms/ml ampicillin at 30°C with shaking overnight. A frozen stock in 40% glycerol was made from this culture and stored at -70°C. DNA was isolated from the culture using Qiagen-tip 100 columns as described above.

[0103] The purified DNA, pTAF172-2, was sequenced using the di-deoxy method and shown to contain the sequence of pLCSF158A coding for M-CSF NA3CA158 with His 9 and His 15 mutated to encode Ala.

[0104] The M-CSF mutein NΔ3CΔ158 H9A, H15A encoded by pTAF172-2 was expressed, purified, refolded to form dimeric protein and assayed essentially as described in U.S. Patent No. 4,929,700 Example 5, using 8M urea as a denaturant and in the DEAE purification step.

[0105] The N-terminal sequence of the purified mutein was determined through 20 cycles, using a standard automated Edman degradation method, and shown to be identical to that of the parental NΔ3CΔ158 M-CSFα reference protein except that His 9 and His 15 had been altered to Ala. Protein concentration was determined using A280 and an extinction co-efficient of 0.68.

[0106] The purified mutein dimers were subjected to bioassay using NFS-60 cells which is an M-CSF dependent cell line which forms colonies in the presence of active M-CSF. Standards and purified mutein samples were serially diluted 1:2 in RPMI media (10% fetal bovine serum) in a 50 microliter volume in 96-well microtiter plates. 50 microliters NFS-60 cells (ATCC NO. CRL 1838), maintained in 4000U/ml M-CSF, washed 2X, and diluted to a concentration of 1 x 10<sup>5</sup> cell/ml, were added to each sample well. Plates were incubated at 37°C, 5% CO<sub>2</sub> for 44h. 5 mg/ml MTT dye was added to each sample and incubation continued for 4h. 20% SDS was added to each sample. The plates were wrapped in foil and left overnight. The absorbance of each sample was read at 570 nm and compared to a standard known M-CSF concentration. The H9A, H15A mutein showed a specific activity 7.6x10<sup>3</sup> U/mg compared to 6.9x10<sup>7</sup> U/mg for the parental M-CSF NΔ3CΔ158 reference in the same assay. This represents a nearly 10,000 fold reduction biological activity for the mutein. The same M-CSF mutein preparation was shown to have greatly decreased M-CSF receptor-binding ability using the NFS-60 receptor competition assay described in this example. Because the H9A, H15A M-CSF mutein was otherwise not significantly different from the parental M-CSFα, including crystallizability and space group (see Example 12) we believe the decrease in biological activity is not due to gross deformation of structure but reflects an alteration in important M-CSF receptor contacts.

[0107] Using essentially the same methodology, two M-CSF muteins contacting singly substituted histidines at residues 9 and 15 were generated (e.g. H9A and H15A). The H9A construct utilized LF80:

5'-CAGGAGAAAGCTTATGTCTGAATATTGTAGCGCCATGATTGGGAGTGGACACCTGCAG-3' (SEQ ID NO 3);

and LF73 (described in this example) as the PCR primers. The H15A construct utilized LF81:

35

5'-CAGGAGAAAGCTTATGTCTGAATATTGTAGCCACATGATTGGGAGTGGAGCCCTGCAG-3' (SEQ ID NO. 4);

and LF 73 (described in this example) as PCR primers. Biological assay of the purified muteins immediately following the refolding step described above showed approximate biological specific activities as follows: 4 x 10<sup>6</sup> U/mg for H9A and less than 3 x 10<sup>3</sup> U/mg for H15A, in an assay which the parental M-CSF construct displayed 8 x 10<sup>7</sup> U/mg. This information, combined with that described above, suggests that H15A as well as possible H9A represents contacts that are important for M-CSF receptor binding. Nearby solvent accessible residues such as Y6 and S13 (see also Table 1) may also represent M-CSF receptor contact residues. Non-proteinaceous mimics of the side chains of H9, H15, and nearby solvent accessible side chains may represent M-CSF agonists or antagonists. Such residues should be left unchanged in M-CSF mutein constructs designed to retain full M-CSF receptor binding but to have M-CSF antagonist properties because they lack significant M-CSF bioactivity. Homodimers of muteins that retain full receptor-binding ability and display significantly reduced bioactivity should represent M-CSF antagonists. M-CSF muteins that are greatly decreased in both M-CSF bioactivity and receptor binding ability (such as H15A) may generally be useful in M-CSF immunoassay applications and might represent useful therapeutic agents for patients having auto-antibodies to M-CSF.

#### Example 11

## Preparation of Q20A, V78K M-CSF Muteins

Using essentially the same methodology described in Example 10, a double mutant of M-CSF (Q20A, V78K) was constructed to test the importance of solvent accessible residues in the central portion of helices A and C. The following PCR primers were used.

#### LF63:

5'-AGGAGAAAGCTTATGTCTGAATATTGTAGCCACATGATTGGGAGTGGACACCTGCAGT CTCTGGCTCGGCTG-3' (SEQ ID NO. 5)

#### LF64:

5'-GGACGCAGGCCTTGTCATGCTCTTCATAATCCTTGGTGAAGCAGCTCTTCAGCCTCAA AGAGAGTTCCTGCAGCTGTTTAATGGC-3' (SEQ ID NO. 6)

20

10

15

[0109] The resulting mutein was expressed, refolded, purified, and assayed as described in Example 8. The specific biological activity was 1.4 X  $10^7$  U/mg, approximately 8-10 fold lower than that of the parental M-CSF $\alpha$  reference standard. The receptor binding activity of this mutein was also decreased.

This result again supports the prediction of over crystallographic study of truncated M-CSFa which con-[0110] cluded that important M-CSF receptor contact residues exist among the solvent accessible residues in helices A and/or C and/or D. Certain of these mutations will, as we have shown, have lower biological activity and lower M-CSF receptorbinding ability. Some may have lower biological activity without a decrease in receptor-binding ability. Some may have increased biological activity and receptor binding ability, and some may have no affect on either.

[0111] Two examples of the latter are Q17A, R21A (produced using PCR primers LF72:

## 5'-TTGTAGCCACATGATTGGGAGTGGACACCTGGCGTCTCTGCAGGCGCTGATTGAC-3'

and LF73 (described in Example 9) and E115A, N119A (produced using LF75:

5 '-CATGACAAGGCCTGCGTCCGAACTTTCTATGAGACACCTCTCCAGTTGCTGGCGAAGG TCAAGGCTGTCTTTAATG-3' (SEQ ID NO. 7):

and

40

45

## LF79: 5'-GGATCAGGATCCCTCGGACTGCCTCTC-3' (SEQ ID NO. 8)).

Both of these constructs changed side chain properties of solvent-accessible amino acids in the areas of interest but did not affect biological specific activity, compared to the parental reference molecule. These results indicate that residues Q17, R21, E115, and N119 do not need to be altered in muteins designed to have M-CSF agonist or antagonist activity. In fact, to minimize the likelihood of antibody formation to potentially administered M-CSF-based proteinaceous drugs, it is desirable to retain the solvent-accessible parental M-CSF residues (to resemble the native molecule) whenever possible.

[0112] The retained activity of the muteins including changes at Q17, R21, E115, and N119 does not rule out large effects on activity contributed by nearby residues (such as H15). In fact, the regions we have altered are predicted by the crystal structure to be important for receptor binding and/or signaling. Antagonistic M-CSF muteins may require use of multiple residue changes per mutein or use of heterodimeric molecules containing one or more mutations in each polypeptide chain, since M-CSF residues important in receptor signaling are believed to be composed of discontinuous regions of M-CSF.

#### Example 12

5

# Formation of M-CSF Heterodimers Having Decreased Receptor-Binding Ability and/or Decreased Biological Specific Activity

[0113] M-CSF can be folded *in vitro* to generate fully active heterodimers, as shown in Example 8. By making heterodimers of M-CSF which incorporate M-CSF muteins with altered M-CSF signaling ability, it should be possible to generate antagonists of M-CSF useful for treatment of patients with M-CSF mediated diseases. To generate a heterodimer containing one subunit of M-CSFαNΔ3/CΔ158, H9A/H15A and one subunit of M-CSFβNΔ3/CΔ221, C157S/C159S, each mutein was expressed in *E. coli* and purified separately by DEAE Sepharose under denaturing and reducing conditions as described in Example 8. The two muteins subunits were mixed together prior to refolding to generate a solution containing a final mutein molar ratio of 1:1, then this solution was diluted to 0.2 mg/ml with refolding buffer as described in Example 8. Following refolding, the heterodimeric molecule was separated from the homodimers by two consecutive passes over a Phenyl TSK-5-PW HPLC column as described in Example 8. No contaminants were detected when the purified heterodimer preparation was examined, by non-reduced SDS-PAGE or size exclusion HPLC using a BioSil SEC250 column (BioRad).

The purified heterodimer was submitted to the NFS60 cell based bioassay described in Example 8. The cal-[0114] culated specific activity was 2.9 x 10<sup>6</sup> U/mg which correlated to a 35-fold reduction as compared to the activity of the parental M-CSF heterodimer described in Example 10. The relative binding affinity to cell surface M-CSF receptor was measured by radioligand displacement in which the displacement of <sup>125</sup>I-M-CSF from an M-CSF receptor by an M-CSF mutein was measured using methods well known in the art. In brief, the following were added in a final volume of 100µl in each well of a 96-well cell culture plate: approximately 80,000 cpm of purified recombinant human M-CSF labeled with <sup>125</sup>I (using lodobeads as described by the manufacturer, Pierce, Rockford, IL), 300,000 NFS-60 cells that had been washed and then grown for 18 hours in growth medium minus the normal maintenance level of M-CSF, plus unlabeled M-CSF that had been serially diluted in the same medium. The plates were incubated at 4°C for 20 hours and the cells were collected on glass-fiber filters. Maximum binding was measured in the absence of unlabeled M-CSF and non-specific binding was measured in the presence of 1000-fold greater concentration of unlabeled M-CSF (compared to labeled M-CSF). The concentration of M-CSF required to inhibit 50% of the labeled M-CSF binding to the cells (IC<sub>50</sub>) was used to determine differences in affinity. Results are expressed as percent displacement of radioactive M-CSF versus mutein concentration (Figure 7). The IC<sub>50</sub> of the heterodimer (Fig. 7 closed squares) was reduced 30-fold to about 500 pM as compared to an IC50 of about 17 pM for M-CSF $\alpha$ N $\Delta$  3/C $\Delta$ 158 (158) (Fig. 7 closed circles). The similarity between the reduction in specific activity and receptor affinity of the heterodimer indicates that the reduction in bioactivity was due to decreased receptor-binding ability. Similarly, the binding affinities of the Q20A, V78KF (Fig. 7 open circles) and H9A/H15A (open squares) muteins were also measured in this radioligand displacement assay. The Q20A,V78K mutein had an IC $_{50}$  of about 100 pM and the H9A,H15A mutein has an IC $_{50}$  of about 1  $\mu$ M, correlating to decreased binding affinities of 5-fold and 50,000-fold, respectively. For each mutein, the reduction in receptor affinity was similar to the reduction in specific activity, again indicating that the reduction in bioactivity was due to reduced receptor-binding ability.

#### 40 Example 13

# Crystallization and Characterization Of M-CSF H9A, H15A Muteins

[0115] The H98, H15A mutein described in Example 10 was crystallized using the hanging drop method described in Examples 1 and 2 using the following buffer conditions: 30% polyethylene glycol 4000; 100 mM Li<sub>2</sub> SO<sub>4</sub>; and 100 mM Tris pH 8.5. The crystals produced under these conditions were rhombohedral prisms having dimensions of 0.7 mm x 0.2 mm x 0.2 mm. X-ray crystallographic analysis using precession photographs showed crystals in the P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub> space group with cell dimensions of a = 33.99, b = 65.37, c = 159.90, d = 90, e = 90, and f = 90 angstroms and diffract to a nominal resolution of 3 angstroms. These physical properties are essentially the same as those observed for the parental N $\Delta$ 3C $\Delta$ 158 M-CSF $\alpha$  molecule and suggests that the biological effects of the H9A, H15A alterations are not the consequence of gross global alterations in M-CSF structure, but rather are the result of altered side chains that are important in interacting with the M-CSF receptor. Alteration of those histidine side chains may have affected receptor binding by changing atoms that interact with, stabilize or facilitate receptor binding or changes in receptor conformation. Changes such as H15A may also have affected these functions by altering the position of the nearby side chain in M-CSF, most likely in the A and/or C helix regions.

[0116] The foregoing examples are presented by way of example and are not intended to limit the scope of the invention as set forth in the appended claims.

# APPENDIX 1

5

55

10	MOTA	10	CA	SER	4	63.753 80.590 222.385 1.00 58.89
	MOTA	13	CA	ŒIJ	5	64.883 77.664 219.972 1.00 59.56
	MOTA	23	CA	TYR	6	62.285 76.840 217.324 1.00 54.54
	MOTA	37	CA	CYS	7	63.509 80.109 215.834 1.00 54.96
15	MOTA	44	CA	SER	8	66.853 79.529 214.160 1.00 54.41
	MOTA	52	CA	HIS	9	65.466 77.798 211.053 1.00 55.11
	MOTA	65	CA	MET	10	61.857 78.767 211.073 1.00 50.66
	MOTA	74	CA	ILE	11	62.173 80.905 207.970 1.00 48.31
20	MOTA	83	CA	GLY	12	63.487 78.618 205.354 1.00 52.83
	MOTA	88	CA	SER	13	64.609 79.967 201.952 1.00 53.13
	ATOM	96	CA	GLY	14	61.665 78.009 200.514 1.00 46.05
05	MOTA	101	CA	HIS	15	59.455 80.924 201.933 1.00 40.25
25	MOTA	114	CA	LEU	16	62.210 83.200 200.992 1.00 41.58
	MOTA	123	CA	GLN	17	62.153 82.266 197.328 1.00 46.17
	MOTA	135	CA	SER	18	58.378 81.868 197.227 1.00 46.55
30	ATOM	143	CA	LEU	19	58.280 85.544 198.199 1.00 44.75
	MOTA	152	CA	GLN	20	60.814 86.044 195.458 1.00 39.32
	MOTA	164	CA	ARG	21	58.457 84.718 192.736 1.00 36.91
	MOTA	181	CA	LEU	22	55.610 86.694 194.252 1.00 38.27
35	MOTA	190	CA	ILE	23	57.676 89.822 193.465 1.00 34.84
	MOTA	199	CA	ASP	24	58.677 88.274 190.086 1.00 31.50
	MOTA	208	CA	SER	25	55.086 87.802 188.978 1.00 31.73
	MOTA	216	CA	GLN	26	54.284 91.540 189.154 1.00 33.24
40	MOTA	228	CA	MET	27	53.961 93.634 186.082 1.00 34.82
	ATOM	237	CA	ŒW	28	56.227 96.566 186.089 1.00 36.54
	ATOM	247	CA	THR	29	53.758 99.303 186.920 1.00 44.80
45	ATOM	256	CA	SER	30	54.590 102.747 188.323 1.00 49.40
	ATOM	264	CA	CYS	31	51.427 103.423 190.046 1.00 43.03
	MOTA	271	CA	GLN	32	52.166 104.336 193.613 1.00 33.74
	MOTA	283	CA	ILE	33	50.430 102.961 196.634 1.00 32.48
50	MOTA	292	CA	THR	34	50.240 104.213 200.133 1.00 38.94
	ATOM	301	CA	PHE	35	51.291 102.173 203.199 1.00 39.16
	MOTA	313	CA	ŒIJ	36	52.707 102.849 206.761 1.00 32.09
	MOTA	323	CA	PHE	37	56.073 101.358 207.299 1.00 25.78

```
MOTA
               335 CA VAL
                               38
                                      59.082 101.751 209.490
                                                                1.00 34.46
       MOTA
               343
                   CA ASP
                               39
                                      61.044 104.798 208.714
                                                                1.00 44.03
 5
       MOTA
               352
                   CA GLN
                               40
                                      64.648 103.690 208.314
                                                                1.00 54.20
       MOTA
               364
                   CA GLU
                               41
                                      65.924 107.142 209.332
                                                                1.00 52.01
       ATOM
               374
                   CA GLN
                               42
                                      63.934 107.629 212.631
                                                                1.00 44.26
       MOTA
               386
                   CA LEU
                                      64.770 104.161 213.955
                               43
                                                                1.00 45.54
 10
               395
       MOTA
                   CA
                       ALA
                               44
                                      68.126 102.952 212.789
                                                                1.00 51.53
       MOTA
               401
                   CA ASP
                               45
                                      69.175 100.197 215.232
                                                                1.00 47.92
       MOTA
               410
                   CA
                       PRO
                               46
                                      68.861 97.054 213.098
                                                               1.00 45.27
15
       MOTA
               417
                   CA
                       VAL
                               47
                                      67.352
                                              94.587 215.613
                                                               1.00 39.92
       MOTA
               425
                   CA
                       CYS
                               48
                                      64.862 .97.220 216.692
                                                               1.00 36.44
       MOTA
               432
                   CA
                       TYR
                               49
                                      64.089
                                              98.158 213.105
                                                               1.00 37.39
       MOTA
               446
                   CA
                       LEU
                               50
                                      63.436
                                              94.474 212.189
                                                               1.00 36.15
20
       MOTA
               455
                   CA LYS
                               51
                                      61.603
                                               94.029 215.472
                                                               1.00 38.66
       MOTA
              468
                   CA LYS
                               52
                                      59.322
                                              96.862 214.367
                                                               1.00 41.72
       MOTA
              481
                   CA ALA
                              53
                                      59.229
                                              96.177 210.534
                                                               1.00 35.81
       MOTA
              487
                       PHE
                   CA
25
                              54
                                      58.096
                                              92.701 211.432
                                                               1.00 40.82
                   CA LEU
      MOTA
              499
                              55
                                      55.236
                                              93.733 213.614
                                                               1.00 42.59
      MOTA
              508
                   CA LEU
                              56
                                      54.205
                                              95.869 210.697
                                                               1.00 44.24
      MOTA
              517
                   CA VAL
                              57
                                      54.596
                                              93.257 207.992
                                                               1.00 35.10
30
      ATOM
              525
                   CA GLN
                              58
                                      51.646
                                              91.431 209.628
                                                               1.00 39.91
      MOTA
              537
                   CA ASP
                              59
                                      49.277
                                              94.254 208.663
                                                               1.00 39.60
      MDTA
              546
                   CA ILE
                              60
                                      50.755
                                              94.729 205.185
                                                               1.00 37.05
      MDTA
              555
                   CA MET
                              61
                                      50.303
                                              90.932 204.676
35
                                                               1.00 33.86
      MDTA
              564
                   CA GLU
                              62
                                      46.766
                                              90.931 206.028
                                                               1.00 48.10
      ATOM
              574
                   CA ASP
                              63
                                      45.602
                                              94.198 204.286
                                                               1.00 49.51
      MOTA
              583
                   CA THR
                              64
                                      47.659
                                              95.244 201.187
                                                               1.00 42.65
40
      MOTA
              592
                   CA MET
                              65
                                      49.106
                                              92.155 199.603
                                                               1.00 37.90
      MOTA
              601
                   CA ARG
                              66
                                      45.646
                                              91.249 198.376
                                                               1.00 39.50
      MOTA
              618
                   CA PHE
                              67
                                      45.516
                                              88.245 195.824
                                                               1.00 35.13
      MOTA
              630
                   CA ARG
                              68
                                      42.655
                                              86.013 194.326
                                                               1.00 44.68
45
      MOTA
              647
                   CA ASP
                              69
                                      41.727
                                              83.454 196.865
                                                               1.00 47.42
      MOTA
              656
                   CA ASN
                              70
                                      43.117
                                              80.057 195.984
                                                               1.00 44.47
      MOTA
              667 CA THR
                              71
                                      45.873 81.088 193.518
                                                               1.00 34.80
```

55

```
MOTA
               676 CA PRO
                              72
                                              80.251 194.114 1.00 29.78
                                      49.607
       MOTA
               683
                    CA ASN
                              73
                                      50.307
                                              83.745 195.418
                                                              1.00 28.47
       MOTA
 5
               694
                   CA ALA
                              74
                                      47.469
                                              83.922 197.953
                                                              1.00 23.24
       MOTA
               700 CA ILE
                              75
                                      48.690
                                              80.431 199.144
                                                               1.00 22.87
       MOTA
               709 CA ALA
                              76
                                              81.892 199.669
                                      52.212
                                                              1.00 29.47
       MOTA
               715
                   CA ILE
                              77
                                      50.816
                                              84.796 201.668
 10
                                                              1.00 37.21
       MOTA
               724
                   CA VAL
                              78
                                      48.902
                                              82.158 203.857
                                                              1.00 33.69
       ATOM
               732
                   CA GLN
                              79
                                      52.207
                                              80.382 204.466
                                                              1.00 33.36
       MOTA
              744
                   CA LEU
                              80
                                     54.040
                                              83.518 205.449
                                                              1.00 32.55
       MOTA
              753
                   CA GIN
 15
                              81
                                     51.152
                                             84.517 207.743
                                                              1.00 35.07
       MOTA
              765
                   CA GLU
                              82
                                     51.541
                                              81.115 209.443
                                                              1.00 37.26
      MOTA
              775
                   CA LEU
                              83
                                     55.367
                                             81.348 209.619
                                                              1.00 35.10
      MOTA
              784
                   CA
                      SER
                              84
                                             84.840 211.066
                                     54.960
                                                              1.00 39.64
      MOTA
20
              792
                   CA LEU
                              85
                                     52.608 83.541 213.704
                                                              1.00 42.55
      MOTA
              108
                   CA
                      ARG
                             86
                                     54.984
                                             80.903 215.121
                                                              1.00 42.96
      ATOM
              818
                   CA
                     LEU
                             87
                                     57.916 83.317 214.670
                                                              1.00 45.11
      MOTA
              827
                   CA
                     LYS
                             88
                                     56.504
                                             85.876 217.138
25
                                                              1.00 50.71
      ATOM
              840
                  CA
                      SER
                             89
                                     57.621
                                             83.610 219.825
                                                              1.00 54.58
      MDTA
              848
                  CA
                      CYS
                             90
                                     61.211
                                             84.524 219.014
                                                              1.00 46.79
      MOTA
              855
                  CA
                     PHE
                             91
                                     60.235
                                             88.070 219.743
                                                              1.00 46.37
      MOTA
              867
                  CA THR
30
                             92
                                     59.619 89.812 223.026
                                                             1.00 52.01
      ATOM
              876
                  CA
                     ALA
                             93
                                     56.916 92.445 223.523
                                                             1.00 56.09
      MOTA
              882
                  CA ASP
                             94
                                     57.914 96.123 223.998
                                                             1.00 59.35
      MOTA
             891
                 CA TYR
                             95
                                    55.685 99.100 225.194
                                                             1.00 66.22
35
      MOTA
             905
                  CA GLU
                             96
                                    52.401 100.066 223.612
                                                             1.00 65.51
      ATOM
             915
                  CA GLU
                             97
                                    53.343 103.890 223.182
                                                             1.00 62.60
     MOTA
             925 CA HIS
                             98
                                    56.046 102.712 220.836
                                                             1.00 58.47
     ATCM
             938
                 CA ASP
                             99
                                    53.422 100.733 218.887
40
                                                             1.00 59.97
     MDTA
             947 CA LYS
                            100
                                    52.950 103.470 216.162
                                                             1.00 57.00
     MOTA
             960
                 CA ALA
                            101
                                    56.259 105.270 216.975
                                                             1.00 48.82
     MOTA
             966
                 CA CYS
                            102
                                    58.861 105.450 214.154
                                                             1.00 42.41
     MOTA
             973 CA VAL
45
                           103
                                    56.175 104.773 211.621
                                                             1.00 33.87
     MOTA
             981
                 CA ARG
                           104
                                    55.993 106.851 208.456
                                                             1.00 42.76
     MOTA
             998 CA THR
                           105
                                    53.577 106.853 205.501
```

50

55

1.00 38.00

```
MOTA
             1007
                   CA PHE
                            106
                                     54.931 105.944 202.136
                                                             1.00 33.48
      MOTA
             1019
                   CA TYR
                            107
                                     53.853 106.546 198.455
                                                              1.00 34.28
 5
      MOTA
             1033
                   CA GLU
                                     55.765 103.906 196.543
                            108
                                                              1.00 33.88
      MOTA
             1043
                  CA THR
                            109
                                     55.266 101.309 193.854
                                                              1.00 36.94
      MOTA
             1052
                  CA PRO
                            110
                                     53.745 97.857 194.657
                                                              1.00 29.84
      MOTA
             1059
                  CA LEU
                            111
                                     57.132 96.374 193.761
 10
                                                              1.00 35.66
      MOTA
             1068
                  CA GLN
                            112
                                     58.991
                                             98.756 196.160
                                                             1.00 33.48
      MOTA
             1080
                  CA LEU
                            113
                                     56.612 97.691 198.912
                                                             1.00 24.26
             1089
      MOTA
                  CA LEU
                            114
                                     57.253 93.932 198.091
                                                             1.00 30.75
 15
      MOTA
             1098
                  CA GLU
                            115
                                     61.063 94.385 197.984
                                                             1.00 39.68
      MOTA
                  CA LYS
             1108
                            116
                                     60.805 95.824 201.526
                                                             1.00 39.24
      MOTA
             1121
                  CA VAL
                            117
                                     58.669 92.840 202.743
                                                             1.00 32.46
      MOTA
             1129
                  CA LYS
                            118
                                     61.312 90.620 201.111
                                                             1.00 37.34
20
      MOTA
             1142
                  CA ASN
                            119
                                     64.231
                                             92.316 202.959
                                                             1.00 39.46
      MOTA
             1153
                  CA VAL
                            120
                                     62.470 92.033 206.370
                                                             1.00 32.07
      MOTA
            1161
                  CA PHE
                            121
                                     61.935
                                             88.268 206.034
                                                             1.00 26.35
      MOTA
             1173 CA ASN
25
                                             88.049 204.416
                            122
                                     65.404
                                                             1.00 36.04
      MOTA
             1184
                 CA GUU
                            123
                                     67.372
                                             89.816 207.138
                                                             1.00 40.82
      MOTA
             1194 CA THR
                            124
                                             88.435 210.122
                                     65.331
                                                             1.00 41.72
      MOTA
             1203 CA LYS
                            125
                                     66.456
                                             85.000 208.648
                                                             1.00 43.07
30
      MOTA
            1216 CA ASN
                            126 -
                                     70.010
                                             86.339 208.179
                                                             1.00 48.89
     MDTA
            1227 CA LEU
                            127
                                     70.000
                                             87.300 211.903
                                                             1.00 47.97
     MOTA
            1236 CA LEU
                            128
                                     68.315
                                             84.194 213.435
                                                             1.00 47.21
     MOTA
            1245
                  CA ASP
35
                            129
                                     71.086 82.134 211.864
                                                             1.00 50.11
     MOTA
            1254
                  CA LYS
                            130
                                     73.448
                                             84.146 214.117
                                                             1.00 52.04
     MOTA
            1267
                  CA ASP
                            131
                                     71.442
                                             84.303 217.451
                                                             1,00 49.58
     MOTA
            1276 CA TRP
                            132
                                     67.946
                                             82.685 218.322
                                                             1.00 53.33
40
     MOTA
            1292
                  CA ASN
                            133
                                     67.568 85.178 221.086
                                                             1.00 59.61
     MOTA
            1303
                  CA ILE
                            134
                                             88.447 219.305
                                     68.572
                                                             1.00 49.82
     MOTA
            1312 CA PHE
                            135
                                     64.945
                                            89.516 219.045
                                                             1.00 54.55
     MDTA
            1324 CA
                      SER
45
                            136
                                     64.745
                                            89.689 222.796
                                                             1.00 62.26
     MOTA
            1332
                  CA LYS
                            137
                                     66.563
                                            93.119 222.569
                                                             1.00 54.38
     MOTA
            1345 CA ASN
                            138
                                    64.785
                                            95.734 224.504
                                                             1.00 58.71
     ATOM
            1356 CA CYS
                            139
                                    63.902 98.337 222.050
                                                             1.00 51.21
```

55

```
MOTA
              1363 CA ASN
                             140
                                     61.704 100.578 224.157 1.00 50.25
       MOTA
              1374 CA ASN
                             141
                                     64.219 103.338 224.061
                                                             1.00 58.30
       MOTA
 5
              1385 CA SER
                             142
                                     65.158 102.992 220.352
                                                             1.00 52.18
       MOTA
              1393 CA PHE
                             143
                                     61.498 103.434 219.645
                                                             1.00 45.87
       MOTA
              1405 CA ALA
                             144
                                     61.829 106.488 221.831
                                                             1.00 47.49
       MOTA
              1411 CA GLU
                             145
                                     64.873 108.022 220.050
                                                             1.00 49.94
 10
       MOTA
              1421 CA CYS
                             146
                                     62.489 108.299 217.017
                                                             1.00 53.65
       MOTA
              1428 CA SER
                             147
                                     60.005 111.197 217.009
                                                             1.00 57.45
       MOTA
              1436 CA SER
                             148
                                     57.163 110.805 214.395
                                                             1.00 59.53
       MOTA
              1444 CA ALA
                             149
                                     54.284 108.658 214.945
                                                             1.00 59.20
 15
       MOTA
              1450 CA GLY
                             150
                                     51.655 108.783 212.248
                                                             1.00 61.95
       ATOM
              1455 CA HIS
                             151
                                     48.764 106.945 210.680
                                                             1.00 69,20
       MOTA
              1468 CA GLU
                             152
                                     46.475 107.377 207.803 1.00 76.61
       MOTA
              1478 CA ALA
                             153
                                     45.813 106.188 204.264
20
                                                             1.00 78.85
       MOTA
              1492 CA SER
                            404
                                     43.875 81.916 155.536
                                                            1.00 38.68
       MOTA
              1495 CA GLU
                            405
                                     41.939
                                            78.947 156.798 1.00 43.57
       MDTA
              1505 CA TYR
                            406
                                     44.549 78.204 159.422
                                                            1.00 41.97
25
       MOTA
             1519 CA CYS
                            407
                                     43.251 81.196 161.434
                                                            1.00 36.99
       MOTA
             1526 CA SER
                            408
                                     41.011 78.917 163.442
                                                            1.00 34.02
       MOTA
             1534 CA HIS
                            409
                                     43.988 77.122 165.021 1.00 40.02
       MOTA
             1547 CA MET
                            410
                                     46.155 79.994 165.895
                                                            1.00 36.48
30
       MOTA
             1556 CA ILE
                            411
                                     44.849 80.843 169.353
                                                            1.00 37.88
       MOTA
             1565 CA GLY
                            412
                                     44.508 77.692 171.141
                                                            1.00 39.09
       MOTA
             1570 CA SER
                            413
                                     43.275 77.480 174.756
                                                            1.00 45.21
       MDTA
             1578 CA GLY
                            414
                                     46.813 77.142 176.120
                                                            1.00 41.00
35
       MOTA
             1583 CA HIS
                            415
                                     47.238 80.801 175.208
                                                            1.00 38.01
       MOTA
             1596 CA LEU
                            416
                                     44.181
                                            81.554 177.458
                                                            1.00 39.12
       MOTA
             1605 CA GLN
                            417
                                     45.501 79.575 180.501
                                                            1.00 39.79
       MOTA
40
             1617 CA SER
                            418
                                     48.547 81.846 180.317 1.00 36.95
      MOTA
             1625 CA LEU
                            419
                                     46.482 85.106 180.212 1.00 33.42
       MOTA
             1634 CA GLN
                            420
                                     44.564 83.708 183.162
                                                            1.00 38.37
      MOTA
             1646 CA ARG
                            421
                                     47.839 82.875 185.041
                                                            1.00 41.95
45
      MOTA
             1663 CA LEU
                            422
                                     49.151
                                            86.514 184.405
                                                            1.00 33.04
      MOTA
             1672 CA ILE
                            423
                                     45.854 87.864 185.738 1.00 25.14
```

50

```
MOTA
              1681 CA ASP
                             424
                                      45.824
                                              85.549 188.861
                                                               1.00 32.93
       MOTA
              1690
                   CA SER
                             425
                                      49.394
                                              86.494 189.904
                                                               1.00 30.92
 5
       MOTA
              1698
                   CA GLN
                             426
                                      48.776 90.276 189.999
                                                               1.00 29.47
       MOTA
              1710
                   CA MET
                             427
                                      48.379
                                              91.899 193.372
                                                              1.00 31.04
       MOTA
              1719
                   CA GW
                             428
                                      44.978
                                              93.418 193.655
                                                              1.00 42.94
       MOTA
              1729
                   CA
                       THR
                             429
 10
                                      45.965
                                              96.999 194.423
                                                              1.00 48.52
       MOTA
              1738
                       SER
                   CA
                             430
                                      43.656 99.720 193.224
                                                              1.00 48.67
       MOTA
              1746
                   CA CYS
                             431
                                      46.604 101.354 191.803
                                                              1.00 52.34
       MOTA
              1753
                   CA
                      ŒИ
                             432
                                      45.813 102.789 188.396
                                                              1.00 50.77
 15
       MOTA
             1765
                   CA
                       ILE
                             433
                                      47.826 103.166 185.221
                                                              1.00 51.17
       MOTA
             1774 CA THR
                             434
                                      46.923 105.057 182.065
                                                              1.00 52.62
       MDTA
             1783
                   CA PHE
                             435
                                      46.641 103.293 178.672
                                                              1.00 50.53
       MOTA
             1795 CA GLU
                             436
                                      45.167 103.968 175.164
                                                              1.00 50.56
20
       MOTA
             1805
                   CA PHE
                             437
                                      42.643 101.271 174.193
                                                              1.00 46.56
       MOTA
             1817 CA VAL
                             438
                                      39.404 101.233 172.116
                                                              1.00 52.14
       MOTA
             1825
                   CA ASP
                             439
                                      36.342 103.107 173.476
                                                              1.00 59.69
      MOTA
             1834
                   CA GLN
25
                             440
                                      33.750 100.386 173.243
                                                              1.00 63.82
      MOTA
             1846
                   CA GLU
                             441
                                      30.764 102.682 173.229
                                                              1.00 64.79
      MOTA
             1856 CA GLN
                             442
                                      31.810 104.569 170.091
                                                              1.00 59.01
      MOTA
             1868
                  CA LEU
                             443
                                     32.994 101.521 168.139
                                                             1.00 55.78
30
      MOTA
             1877
                   CA LYS
                             444
                                     30.279 99.137 169.086
                                                              1.00 54.82
      MOTA
             1890
                   CA ASP
                             445
                                     30.862 96.668 166.093
                                                              1.00 51.31
      MOTA
             1899
                   CA PRO .
                             446
                                     32.144
                                             93.267 167.135
                                                              1.00 47.34
      MOTA
             1906
                   CA VAL
                             447
35
                                     34.363 92.592 164.111
                                                              1.00 41.42
      MOTA
                   CA CYS
             1914
                             448
                                     35.680 96.184 163.845
                                                              1.00 39.44
      MOTA
             1921
                   CA TYR
                             449
                                     36.273
                                             96.949 167.607
                                                              1.00 40.33
      MOTA
             1935
                   CA LEU
                             450
                                     38.224
                                             93.580 167.624
                                                              1.00 30.87
40
      MOTA
             1944
                   CA LYS
                             451
                                     39.989
                                             94.664 164.462
                                                              1.00 34.49
      MOTA
             1957
                   CA LYS
                             452
                                     40.962
                                             97.918 166.163
                                                             1.00 35.38
      MOTA
             1970 CA ALA
                             453
                                     41.503
                                             96.317 169.572
                                                             1.00 28.00
      MOTA
             1976 CA PHE
                             454
                                     44.012
                                             93.845 168.164
                                                              1.00 32.56
45
      MOTA
             1988 CA LEU
                             455
                                     46.236
                                             96.653 166.861
                                                              1.00 40.19
      MOTA
             1997
                   CA LEU
                            456
                                     46.402
                                             98.605 170.114
                                                             1.00 42.45
      MOTA
             2006 CA VAL
                            457
                                     46.966
                                             95.358 172.060 1.00 38.09
```

55

```
MOTA
              2014 CA GEN
                             458
                                      50.298
                                             95.232 170.161
                                                              1.00 48.02
       MOTA
              2026 CA ASP
                             459
                                      51.630
                                             98.529 171.587
                                                              1.00 41.96
 5
       MOTA
              2035 CA ILE
                             460
                                      50.230
                                             97.686 175.001
                                                              1.00 25.82
       MOTA
              2044 CA MET
                             461
                                      52.409
                                             94.542 175.085
                                                              1.00 33.69
       MOTA
              2053 CA GU
                             462
                                      55.702
                                             96.253 173.969
                                                              1.00 45.36
       MOTA
              2063
                   CA ASP
                             463
                                      55.069
                                             99.093 176.249
                                                              1.00 48.41
 10
       MOTA
              2072 CA THR
                             464
                                             98.555 179.186
                                     52.771
                                                             1.00 41.63
       MOTA
              2081 CA MET
                             465
                                     52.743
                                             94.854 180.171
                                                              1.00 42.05
       MOTA
              2090 CA ARG
                             466
                                     56.333
                                             94.928 181.338
                                                             1.00 44.99
       MOTA
              2107 CA PHE
                             467
                                     57.618
                                             91.923 183.501
                                                             1.00 38.72
 15
       MOTA
              2119 CA ARG
                             468
                                     61.258
                                             90.894 184.513
                                                             1.00 50.64
       MOTA
              2136 CA ASP
                             469
                                     63.046
                                             88.868 181.895
                                                             1.00 53.46
       MOTA
              2145 CA ASN
                             470
                                     62.896 85.096 181.832
                                                             1.00 52.84
       MOTA
              2156 CA THR
20
                             471
                                     60.211
                                             84.798 184.479
                                                             1.00 47.58
       MOTA
              2165 CA PRO
                             472
                                     57.055
                                             82.652 183.744
                                                             1.00 46.82
       MOTA
              2172 CA ASN
                             473
                                     54.971
                                             85.777 183.172
                                                             1.00 38.75
       MOTA
             2183 CA ALA
                             474
                                     57.431
                                             87.402 180.828
                                                             1.00 37.27
25
       MOTA
             2189
                   CA
                      ILE
                             475
                                     57.795
                                             84.219 178.807
                                                             1.00 32.85
             2198 CA ALA
       MOTA
                             476
                                     54.021
                                             84.034 178.455
                                                             1.00 27.98
       MDTA
             2204
                   CA ILE
                             477
                                     53.844
                                             87.522 176.756
                                                             1.00 30.47
       MOTA
             2213
                   CA VAL
                             478
                                     56.690
                                             86.495 174.201
                                                             1.00 30.40
30
       MOTA
             2221 CA GLN
                             479
                                     54.381
                                             83.711 173.325
                                                             1.00 31.18
       MOTA
             2233
                   CA LEU
                             480
                                     51.419
                                             86.038 172.893
                                                             1.00 27.44
       MOTA
             2242 CA GLN
                             481
                                     53.392
                                             88.551 170.821
                                                             1.00 31.59
       MOTA
             2254
                   CA GIU
                             482
                                             85.683 168.506
35
                                     54.572
                                                             1.00 29.90
       MOTA
             2264
                   CA LEU
                             483
                                     50.900
                                             84.532 168.340
                                                             1.00 22.16
       MDTA
             2273
                   CA SER
                             484
                                     50.077
                                             88.156 167.451
                                                             1.00 29.37
       MOTA
             2281 CA LEU
                             485
                                             87.880 164.490
                                     52.410
                                                             1.00 32.70
40
       MOTA
             2290 CA ARG
                             486
                                     50.865
                                             84.572 163.229
                                                             1.00 30.59
       MOTA
             2307
                   CA LEU
                             487
                                     47.396
                                             86.383 163.343
                                                             1.00 31.02
       MOTA
             2316
                  CA LYS
                             488
                                     48.292
                                             89.399 161.352
                                                             1.00 34.90
       MOTA
             2329
                   CA SER
                             489
                                             87.397 158.200
                                     47.728
                                                             1.00 35.68
45
       MOTA
             2337
                   CA
                       CYS
                             490
                                     44.112
                                             87.134 159.417
                                                             1.00 35.03
       MOTA
             2344 CA PHE
                             491
```

50

55

43.520 90.786 159.602

1.00 42.73

```
MOTA
              2356 CA THR
                              492
                                      43.145
                                              93.174 156.663
                                                              1.00 56.24
       MOTA
              2365 CA ALA
                              493
                                      44.877
                                              96.467 156.617
                                                              1.00 57.56
 5
       MOTA
              2371
                    CA ASP
                              494
                                      42.318 99.142 157.302
                                                               1.00 57.16
       MOTA
              2380
                   CA ALA
                             495
                                      42.063 102.909 156.516
                                                              1.00 62.14
       MOTA
              2386 CA GLU
                             496
                                      45.542 104.198 157.805
                                                              1.00 64.30
       MOTA
              2396 CA GLU
                                      44.589 107.547 159.382
                             497
 10
                                                              1.00 70.67
       MOTA
              2406 CA HIS
                             498
                                      41.577 105.660 160.752
                                                              1.00 75.16
       MOTA
              2419 CA ASP
                             499
                                      44.400 103.329 161.968
                                                              1.00 72.11
       MOTA
              2428
                   CA LYS
                             500
                                      44.773 105.712 164.874
                                                              1.00 65.63
 15
       MOTA
              2441
                   CA ALA
                             501
                                      41.040 106.628 165.617
                                                              1.00 55.63
       MOTA
              2447
                   CA CYS
                             502
                                      38.563 105.427 168.336
                                                              1.00 53.33
       MOTA
             2454
                   CA VAL
                             503
                                      41.155 105.639 171.100
                                                              1.00 52.46
       MOTA
              2462 CA ARG
                             504
                                      39.866 106.944 174.425
20
                                                              1.00 51.00
       MOTA
             2479 CA THR
                             505
                                      42.588 106.921 177.117
                                                              1.00 50.42
      -ATOM
             2488 CA PHE
                             506
                                     41.955 104.935 180.331
                                                              1.00 43.64
       MDTA
             2500
                   CA TYR
                             507
                                      43.306 105.190 183.809
                                                              1.00 41.84
       MOTA
             2514
25
                   CA GLU
                             508
                                     42.651 101.611 185.218
                                                              1.00 39.35
       MDTA
             2524
                   CA THR
                             509
                                     44.156 99.096 187.597
                                                              1.00 41.67
       MOTA
             2533 CA PRO
                             510
                                     46.618 96.563 185.969
                                                              1.00 39.70
       MOTA
             2540 CA LEU
                             511
                                     44.134
                                             93.600 186.416
                                                              1.00 42.72
30
      MDTA
             2549 CA GLN
                             512
                                     41.460 95.632 184.624
                                                              1.00 37.47
      MOTA
             2561
                   CA LEU
                             513
                                     43.895 96.148 181.765
                                                              1.00 27.02
      MOTA
             2570 CA LEU
                             514
                                     44.769 92.421 181.627
                                                             1.00 24.55
      ATOM
             2579 CA GLU
                             515
                                     41.076 91.327 181.688
35
                                                             1.00 31.93
      MOTA
             2589 CA LYS
                             516
                                     40.903
                                             93.442 178.515
                                                             1.00 31.96
      ATOM
             2602 CA VAL
                             517
                                     43.981 91.949 176.692
                                                             1.00 26.46
      MOTA
             2610 CA LYS
                             518
                                     42.528
                                             88.538 177.852
                                                             1.00 30.83
40
      ATOM
             2623 · CA ASN
                            519
                                     39.106 89.397 176.432
                                                             1.00 37.62
      MOTA
             2634 CA VAL
                            520
                                     40.221 . 90.488 172.933
                                                             1.00 35.75
      MOTA
             2642 CA PHE
                            521
                                     42.191 87.109 172.662
                                                             1.00 25.45
      MOTA
             2654
                   CA ASN
                            522
                                     39.292 85.079 173.953
45
                                                             1.00 26.73
      MOTA
             2665
                   CA GLU
                            523
                                     36.870
                                             86.722 171.437
                                                             1.00 33.04
      MOTA
             2675
                   CA
                      THR
                            524
                                             86.923 168.369
                                     39.260
                                                             1.00 29.70
      ATOM
             2684 CA LYS
                            525
                                     39.423
                                             83.195 169.011
                                                             1.00 28.41
```

55

```
MOTA
             2697 CA ASN
                             526
                                     35.624
                                             82.663 169.124 1.00 34.46
       MOTA
             2708 CA LEU
                             527
                                     34.816
                                             84.653 166.021
                                                              1.00 28.09
5
       MDTA
             2717
                   CA LEU
                             528
                                     37.617
                                             82.902 164.068
                                                              1.00 25.62
       MOTA
             2726 CA ASP
                             529
                                     35.899
                                             79.662 164.964
                                                              1.00 35.00
       MOTA
             2735 CA LYS
                             530
                                     32.659
                                             81.068 163.537
                                                              1.00 37.54
       MOTA
             2748 CA ASP
                             531
                                     34.250
                                             82.150 160.272
10
                                                              1.00 37.92
       MOTA
             2757 CA TRP
                             532
                                     38.053
                                             82.453 159.573
                                                              1.00 40.02
       MOTA
             2773
                   CA ASN
                             533
                                     38.395
                                             85.446 157.159
                                                              1.00 42.64
       MOTA
             2784
                   CA ILE
                             534
                                     35.840
                                             87.437 159.393
                                                             1.00 39.27
       MOTA
15
             2793
                   CA PHE
                             535
                                     38.635
                                             89.996 159.881
                                                              1.00 40.49
      MOTA
             2805 CA SER
                             536
                                     38.485
                                             90.720 156.238
                                                             1.00 46.56
      MOTA
             2813 CA LYS
                             537
                                     35.844 93.371 157.102
                                                              1.00 45.61
      MOTA
             2826 CA 'ASN
                             538
                                     35.653 96.936 155.760
                                                             1.00 54.39
20
      MOTA
             2837 CA CYS
                             539
                                     35.786 98.626 159.004
                                                             1.00 50.55
      ATOM
             2844 CA ASN
                             540
                                     36.738 101.888 157.286
                                                             1.00 58.49
      MOTA
             2855
                   CA ALA
                             541
                                     33.143 102.934 157.648
                                                             1.00 61.90
      MOTA
             2861
                   CA SER
                            542
                                     32.854 101.855 161.326
                                                             1.00 56.62
25
      MOTA
             2869
                   CA PHE
                             543
                                     35.911 103.888 162.235
                                                             1.00 55.25
      MOTA
             2881
                   CA ALA
                            544
                                     33.885 107.102 161.741
                                                             1.00 61.94
      MOTA
             2887
                   CA GLU
                            545
                                     30.949 107.003 164.203
                                                             1.00 63.83
      MOTA
             2897
30
                   CA CYS
                             546
                                     33.917 107.418 166.572
                                                             1.00 61.90
      MOTA
             2904 CA SER
                            547
                                     33.751 110.878 168.150
                                                             1.00 66.22
      MDTA
             2912 CA ALA
                            548
                                     37.461 110.300 168.811
                                                             1.00 65.91
      END
```

55

50

35

40

## SEQUENCE LISTING

5	(1) GENERAL INFORMATION:
·	(i) APPLICANTS: CHIRON CORPORATION AND THE REGENTS OF THE UNIVERSITY OF CALIFORNIA
	(ii) TITLE OF INVENTION: CRYSTALLIZATION OF M-CSF
10	(iii) NUMBER OF SEQUENCES: 10
15	<ul> <li>(iv) CORRESPONDENCE ADDRESS:</li> <li>(A) ADDRESSEE: Marshall, O'Toole, Gerstein, Murray &amp; Boru</li> <li>(B) STREET: 6300 Sears Tower, 233 South Wacker Drive</li> <li>(C) CITY: Chicago</li> <li>(D) STATE: Illinois</li> </ul>
	(B) COUNTRY: USA (F) ZIP: 60606-6402
20	(v) COMPUTER READABLE FORM:  (A) MEDIUM TYPE: Floppy disk  (B) COMPUTER: IBM PC compatible  (C) OPERATING SYSTEM: PC-DOS/MS-DOS  (D) SOFTWARE: PatentIn Release #1.0, Version #1.25
25	(VI) CURRENT APPLICATION DATA:  (A) APPLICATION NUMBER:  (B) FILING DATE:  (C) CLASSIFICATION:
	<ul><li>(vii) PRIOR APPLICATION DATA:</li><li>(A) APPLICATION NUMBER: US 07/896,512 CIP</li><li>(B) FILING DATE: 09-JUN-1992</li></ul>
30	(viii) ATTORNEY/AGENT INFORMATION:  (A) NAME: Clough, David W.  (B) REGISTRATION NUMBER: 36,107  (C) REFERENCE/DOCKET NUMBER: 31436
35	(ix) TELECOMMUNICATION INFORMATION:  (A) TELEPHONE: 312/474-6300  (B) TELEFAX: 312/474-0448  (C) TELEX: 25-3856
40	(2) INFORMATION FOR SEQ ID NO:1:  (1) SEQUENCE CHARACTERISTICS:
45	(A) LENGTH: 56 base pairs (B) TYPE: nucleic acid (C) STRANDEDNESS: single (D) TOPOLOGY: linear
40	(ii) MOLECULE TYPE: DNA
50	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:1:  AGGTGTCTCA TAGAAAGTTC GGACGCAGGC CTTGTCATGC TCTTCATAAT CCTTGG 50
	(2) INFORMATION FOR SEQ ID NO:2:

29

5	(1) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 58 bas pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
	(ii) MOLECULE TYPE: DNA	
10	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:2:	
	CAGGAGAAAG CTTATGTCTG AATATTGTAG CGCCATGATT GGGAGTGGAG CCCTGCAG	58
	(2) INFORMATION FOR SEQ ID NO:3:	
15	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 58 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
20	(ii) MOLECULE TYPE: DNA	
	(x1) SEQUENCE DESCRIPTION: SEQ ID NO:3:	
25	CAGGAGAAAG CTTATGTCTG AATATTGTAG CGCCATGATT GGGAGTGGAC ACCTGCAG	58
	(2) INFORMATION FOR SEQ ID NO:4:	
30	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 58 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
	(ii) MOLECULE TYPE: DNA	
35	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:4:	
	CAGGAGAAAG CTTATGTCTG AATATTGTAG CCACATGATT GGGAGTGGAG CCCTGCAG	58
40	(2) INFORMATION FOR SEQ ID NO:5:	
	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 72 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
45	(ii) MOLECULE TYPE: DNA	
50	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:5:	
	AGGAGAAAGC TTATGTCTGA ATATTGTAGC CACATGATTG GGAGTGGACA CCTGCAGTCT	60
	CTGGCTCGGC TG	72

(2) INFORMATION POR SEQ ID NO:6:

5	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 85 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
10	(ii) MOLECULE TYPE: DNA	
	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:6:	
	GGACGCAGGC CTTGTCATGC TCTTCATAAT CCTTGGTGAA GCAGCTCTTC AGCCTCAAAG	60
15	AGAGTTCCTG CAGCTGTTTA ATGGC 85	
	(2) INFORMATION FOR SEQ ID NO:7:	
20	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 76 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
25	(ii) MOLECULE TYPE: DNA	
	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:7:	
	CATGACAAGG CCTGCGTCCG AACTTTCTAT GAGACACCTC TCCAGTTGCT GGCGAAGGTC	60
30	AAGGCTGTCT TTAATG 76	
	(2) INFORMATION FOR SEQ ID NO:8:	
35	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 27 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
40	(ii) MOLECULE TYPE: DNA	
	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:8:	
<b>4</b> 5	GGATCAGGAT CCCTCGGACT GCCTCTC	. 27
	(2) INFORMATION FOR SEQ ID NO:9:	
50	(i) SEQUENCE CHARACTERISTICS:  (A) LENGTH: 26 base pairs  (B) TYPE: nucleic acid  (C) STRANDEDNESS: single  (D) TOPOLOGY: linear	
	(ii) MOLECULE TYPE: DNA	

31

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:9:

#### GCGTACCATG GGCCCAGGAG TTCTGC

26

- (2) INFORMATION FOR SEQ ID NO:10:
  - (i) SEQUENCE CHARACTERISTICS:
    - (A) LENGTH: 33 base pairs
    - (B) TYPE: nucleic acid
    - (C) STRANDEDNESS: single
    - (D) TOPOLOGY: linear
  - (ii) MOLECULE TYPE: DNA

15

5

10

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:10:

. AGTCGAGGAT CCTCAATCCG GGGGATGCGT GTG

33

20

#### Claims

- A composition comprising crystalline M-CSF, the M-CSF being a dimer of two M-CSF polypeptide monomers, wherein each monomer has substantially the same amino acid sequence as a polypeptide selected from the group consisting of mature M-CSFα polypeptide, mature M-CSFβ polypeptide, and mature M-CSFγ polypeptide.
  - 2. The composition of claim 1 wherein the crystalline dimer has P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub> symmetry.

3. The composition of claim 2 wherein the M-CSF is biologically active when resolubilized.

- 4. The composition of claim 1 wherein the two M-CSF polypeptide monomers are disulfide-linked to each other.
- 5. The composition of claim 1 wherein each of the two M-CSF polypeptide monomers is an M-CSF $\alpha$  polypeptide.
  - 6. The composition of claim 1 wherein each monomer has an amino acid sequence between about 145 and about 156 amino acids in length, and wherein said sequence has a carboxylic acid terminus located between about 145 and 156 residues measured from the amino terminus of a mature M-CSF polypeptide.

40

30

- 7. The composition of claim 6 wherein at least one of the polypeptides contains residues 4 to 158 of mature M-CSF $\alpha$  polypeptide.
- 8. The composition of claim 1 wherein at least one of the M-CSF polypeptide monomers is unglycosylated.

45

- A method of crystallizing a M-CSFα dimer having two polypeptide monomers, each containing 146 to 162 amino acids of mature M-CSFα, the method comprising the steps of:
  - a) mixing a solution of the M-CSF $\alpha$  dimer with a precipitant, whereby an M-CSF $\alpha$  mixture is formed;
  - b) precipitating crystalline M-CSF $\alpha$  dimer; and
  - c) isolating the crystalline M-CSF $\alpha$  dimer.
- 10. The method of claim 9 wherein the M-CSF $\alpha$  dimer is produced from recombinant M-CSF $\alpha$  isolated from a bacterial cell.

55

- 11. The method of claim 9 wherein the precipitant comprises polyethylene glycol.
- 12. The method of claim 9 wherein at least one of the polypeptide monomers is unglycosylated.

- 13. The method of claim 9 wherein at least one of the polypeptide monomers contains residues 4 to 158 of mature M-CSF $\alpha$  polypeptide.
- 14. The method of claim 9 wherein the crystalline M-CSF $\alpha$  dimer is biologically active when resolubilized.
- 15. The method of claim 9 wherein the step of precipitating the M-CSF $\alpha$  crystal comprises equilibrating the M-CSF $\alpha$  dimer mixture with a precipitant solution.
- 16. The method of claim 15 wherein the precipitant solution comprises a higher concentration of precipitant than the M-CSFα mixture.
  - 17. The method of claim 15 wherein the step of equilibrating comprises applying the M-CSFα mixture to a surface and equilibrating the applied M-CSFα mixture with a reservoir of the precipitant solution.
- 18. The method of claim 9 wherein the step of precipitating the crystalline M-CSFα dimer comprises a step of applying M-CSFα seed crystals to the M-CSFα dimer mixture.
  - 19. The method of claim 9 wherein the step of precipitating the M-CSF $\alpha$  dimer comprises the step of altering the temperature of the M-CSF $\alpha$  mixture.
  - 20. Crystalline M-CSFα (4-158)

5

20

25

35

40

45

50

- 21. The crystalline M-CSF of claim 20 comprising crystals of M-CSF having measurements of at least about 0.3 mm in at least two dimensions.
- 22. The crystalline M-CSF of claim 20 which is capable of diffracting X-ray radiation to produce a diffraction pattern representing the three-dimensional structure of the M-CSF.
- 23. The crystalline M-CSF of claim 22 wherein the diffraction pattern represents the three-dimensional structure of a portion of the M-CSFα (4-158) to an accuracy of at least five angstroms, wherein the portion includes amino acid residues between residue 4 and about residue 153.
  - 24. A method for determining a three-dimensional structure of an M-CSF dimer containing two M-CSF polypeptide monomers, each having between about 146 and 162 amino acids beginning from about residue 1 to about residue 5 of mature M-CSF, the method comprising:
    - a) crystallizing the M-CSF dimer;
    - b) irradiating the crystalline M-CSF to obtain a diffraction pattern characteristic of the crystalline M-CSF; and
    - c) transforming the diffraction pattern into the three-dimensional structure of the M-CSF dimer.
  - 25. A method of using a three-dimensional structure of M-CSF derived from an M-CSF crystal wherein the three-dimensional structure of M-CSF includes an M-CSF receptor-binding region, the method comprising identifying compounds having structures that interact with a receptor-binding region of the three-dimensional structure of M-CSF and function as an M-CSF agonist or antagonist.
  - 26. The method of claim 25 wherein the three-dimensional structure of M-CSF includes alpha-carbon coordinates substantially the same as those of the structural information presented in Appendix 1.
  - 27. A method of identifying M-CSF agonists or antagonists, the method comprising the following steps:
    - a) crystallizing an M-CSF dimer to form at least one M-CSF crystal, the M-CSF dimer containing a group of amino acid residues defining an M-CSF receptor-binding region:
    - b) irradiating an M-CSF crystal produced by the procedure of step (a) to obtain a diffraction pattern of the M-CSF crystal;
    - c) determining the three-dimensional structure of M-CSF from the diffraction pattern, the three-dimensional structure including an M-CSF receptor-binding region; and
      - d) identifying an M-CSF agonist or antagonist compound having a three-dimensional structure that functionally duplicates essential M-CSF receptor-binding solvent accessible residues present in the three-dimensional

structure of the M-CSF receptor-binding region, said M-CSF agonist or antagonist having altered signal transdiction capacity to M-CSF responsive cells.

28. The method of claim 27 wherein said solvent accessible residues do not participate in formation of the M-CSF dimer interface.

#### Patentansprüche

5

15

25

35

40

45

50

- Zusammensetzung, umfassend kristallines M-CSF, wobei M-CSF ein Dimer aus zwei M-CSF-Polypeptidmonomeren ist, in dem jedes Monomer im Wesentlichen die gleiche Aminosäuresequenz wie ein Polypeptid hat, das aus der Gruppe, bestehend aus reifem M-CSFα-Polypeptid, reifem M-CSFβ-Polypetid und reifem M-CSFγ-Polypeptid, ausgewählt wird.
  - 2. Zusammensetzung nach Anspruch 1, worin das kristalline Dimer P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub>-Symmetrie aufweist.
  - 3. Zusammensetzung nach Anspruch 2, worin M-CSF biologisch aktiv ist, wenn es wieder in Lösung gebracht wird.
  - 4. Zusammensetzung nach Anspruch 1, worin die zwei M-CSF-Polypeptidmonomeren disulfidverbrückt sind.
- Zusammensetzung nach Anspruch 1, worin jedes der zwei M-CSF-Polypeptidmonomeren ein M-CSFα-Polypeptid
  ist.
  - 6. Zusammensetzung nach Anspruch 1, worin jedes Monomer eine Aminosäuresequenz einer Länge zwischen etwa 145 und etwa 156 Aminosäuren aufweist und wobei die Sequenz einen Carboxylterminus aufweist, der sich etwa zwischen den Resten 145 und 156 befindet, gemessen von dem Aminoterminus des reifen M-CSF-Polypeptids.
    - Zusammensetzung nach Anspruch 6, worin mindestens eines der Polypeptide die Reste 4 bis 158 des reifen M-CSFα-Polypeptids enthält.
- Zusammensetzung nach Anspruch 1, wobei mindestens eines der M-CSF-Polypeptidmonomeren nicht glykosyliert ist.
  - 9. Verfahren zum Kristallisieren eines M-CSFα-Dimers mit zwei Polypeptidmonomeren, wobei jedes 146 bis 162 Aminosäuren des reifen M-CSFα enthält und das Verfahren die folgenden Schritte umfasst:
    - a) Mischen einer Lösung des M-CSF $\alpha$ -Dimers mit einem Fällungsmittel, wobei ein M-CSF $\alpha$ -Gemisch gebildet wird,
    - b) Ausfällen des kristallinen M-CSF $\alpha$ -Dimers und
    - c) Isolieren des kristallinen M-CSFα-Dimers.
  - Verfahren nach Anspruch 9, worin das Fällungsmittel Polyethylenglycol umfasst.
  - 11. Verfahren nach Anspruch 9, worin das M-CSF $\alpha$ -Dimer aus rekombinantem, aus einer Bakterienzelle isoliertem M-CSF $\alpha$  hergestellt wird.
  - 12. Verfahren nach Anspruch 9, worin mindestens eines der Polypeptidmonomere nicht glykosyliert ist.
  - 13. Verfahren nach Anspruch 9, worin mindestens eines der Polypeptidmonomere die Reste 4 bis 158 des reifen M- $CSF\alpha$ -Polypeptids enthält.
  - 14. Verfahren nach Anspruch 9, worin das kristalline M-CSF $\alpha$ -Dimer biologisch aktiv ist; wenn es wieder in Lösung gebracht wird.
  - Verfahren nach Anspruch 9, worin der Fällungsschritt des M-CSFα-Kristalls das Equilibrieren des M-CSFα-Dimergemisches mit einer Fällungsmittel-Lösung umfasst.
  - Verfahren nach Anspruch 15, in dem die Lösung des Fällungsmittels eine h\u00f6here F\u00e4llungsmittelkonzentration umfasst als das M-CSFα-Gemisch.

- 17. Verfahren nach Anspruch 15, in dem der Equilibrierungsschritt das Auftragen des M-CSFα-Gemisches auf eine Oberfläche und das Equilibrieren des aufgetragenen M-CSFα-Gemisches mit einem Reservoir der Lösung des Fällungsmittels umfasst.
- Verfahren nach Anspruch 9, in dem der Fällungsschritt des kristallinen M-CSFα-Dimers einen Schritt umfasst, bei dem M-CSFα-Impfkristalle dem M-CSFα-Dimergemisch zugeführt werden.
  - Verfahren nach Anspruch 9, in dem der Fällungsschritt des M-CSFα-Dimers einen Schritt umfasst, bei dem die Temperatur des M-CSFα-Gemisches verändert wird.
  - 20. Kristallines M-CSFa (4-158).

10

15

25

30

35

45

- Kristallines M-CSF nach Anspruch 20, umfassend M-CSF-Kristalle mit Abmessungen von mindestens etwa 0,3
  mm in mindestens zwei Dimensionen.
- 22. Kristallines M-CSF nach Anspruch 20, das in der Lage ist, Röntgenstrahlung zu beugen, um ein Beugungsbild zu erzeugen, das die dreidimensionale Struktur von M-CSF repräsentiert.
- 23. Kristallines M-CSF nach Anspruch 22, wobei das Beugungsbild die dreidimensionale Struktur eines Teils von M-CSFα (4-158) mit einer Genauigkeit von mindestens 5 Angström repräsentiert, und dieser Teil Aminosäurereste zwischen Rest 4 und etwa Rest 153 umfasst.
  - 24. Verfahren zur Bestimmung einer dreidimensionalen Struktur eines M-CSF-Dimers mit zwei M-CSF-Polypeptidmonomeren, wobei jedes etwa zwischen 146 und 162 Aminosäuren aufweist, ausgehend von etwa Rest 1 bis etwa Rest 5 des reifen M-CSF, das Verfahren umfassend:
    - a) Kristallisieren des M-CSF-Dimers;
    - b) Bestrahlen von kristallinem M-CSF, um ein für kristallines M-CSF charakteristisches Beugungsbild zu erhalten und
    - c) Transformieren des Beugungsbildes in eine dreidimensionale Struktur des M-CSF-Dimers.
  - 25. Verfahren zur Verwendung einer aus einem M-CSF-Kristall abgeleiteten dreidimensionalen Struktur von M-CSF, wobei die dreidimensionale Struktur von M-CSF einen M-CSF-Rezeptorbindungsbereich umfasst und das Verfahren die I-dentifizierung von Verbindungen umfasst, die Strukturen haben, die mit einem Rezeptorbindungsbereich der dreidimensionalen Struktur von M-CSF wechselwirken und als M-CSF-Agonist oder -Antagonist fungieren.
  - 26. Verfahren nach Anspruch 25, wobei die dreidimensionale Struktur von M-CSF alpha-Kohlenstoffkoordinaten umfasst, die im Wesentlichen die gleichen sind wie diejenigen der Strukturinformationen in Anhang 1.
- 40 27. Verfahren zur Identifizierung von M-CSF-Agonisten oder -Antagonisten, wobei das Verfahren die folgenden Schritte umfasst:
  - a) Kristallisieren eines M-CSF-Dimers unter Bildung von mindestens einem M-CSF-Kristall, wobei das M-CSF-Dimer eine Gruppe von Aminosäureresten enthält, die einen M-CSF-Rezeptorbindungsbereich definieren;
  - b) Bestrahlen eines durch das Verfahren von Schritt (a) hergestellten M-CSF-Kristalls, um ein Beugungsbild des M-CSF-Kristalls zu erhalten;
  - c) Bestimmung der dreidimensionalen Struktur von M-CSF aus dem Beugungsbild, wobei die dreidimensionale Struktur einen M-CSF-Rezeptorbindungsbereich umfasst und
  - d) Identifizierung einer M-CSF-Agonisten- oder Antagonistenverbindung mit einer dreidimensionalen Struktur, die in funktioneller Weise wesentliche, dem Lösungsmittel zugängliche M-CSF-Rezeptorbindungsreste, die in der dreidimensionalen Struktur des M-CSF-Rezeptorbindungsbereichs vorliegen, dupliziert, wobei der M-CSF-Agonist oder -Antagonist ein verändertes Signalübertragungsvermögen gegenüber auf M-CSF ansprechenden Zellen aufweist.
- 28. Verfahren nach Anspruch 27, wobei die dem Lösungsmittel zugänglichen Reste nicht an der Bildung der M-CSF-Dimergrenzfläche teilnehmen.

#### Revendications

5

20

30

35

40

- Composition comprenant du M-CSF cristallin, le M-CSF étant un dimère de deux monomères polypeptidiques de M-CSF, où chaque monomère a sensiblement la même séquence d'acides aminés qu'un polypeptide choisi dans le groupe consistant en le polypeptide M-CSFα mature, le polypeptide M-CSFβ mature et le polypeptide M-CSFγ mature.
- 2. Composition selon la revendication 1 où le dimère cristallin a la symétrie P2<sub>1</sub>2<sub>1</sub>2<sub>1</sub>.
- Composition selon la revendication 2 où le M-CSF est biologiquement actif quand il est resolubilisé.
  - Composition selon la revendication 1 où les deux monomères polypeptidiques de M-CSF sont liés l'un à l'autre par une liaison disulfure.
- Composition selon la revendication 1 où chacun des deux monomères polypeptidiques de M-CSF est un polypeptide M-CSFα.
  - 6. Composition selon la revendication 1 où chaque monomère a une séquence d'acides aminés d'une longueur située entre environ 145 et environ 156 acides aminés, et où ladite séquence a une extrémité carboxy-terminale située à environ 145 à 156 résidus en mesurant à partir de l'extrémité amino-terminale d'un polypeptide M-CSF mature.
  - 7. Composition selon la revendication 6 où au moins l'un des polypeptides contient les résidus 4 à 158 du polypeptide M-CSF $\alpha$  mature.
- 25 8. Composition selon la revendication 1 où au moins l'un des monomères polypeptidiques de M-CSF est non glyco-sylé.
  - Procédé de cristallisation d'un dimère de M-CSFα ayant deux monomères polypeptidiques contenant chacun 146 à 162 acides aminés du M-CSFα mature, le procédé comprenant les étapes de:
    - a) mélange d'une solution du dimère de M-CSF $\alpha$  avec un précipitant, de sorte qu'il se forme un mélange de M-CSF $\alpha$ ;
    - b) précipitation d'un dimère de M-CSF $\alpha$  cristallin; et
    - c) isolement du dimère de M-CSF $\alpha$  cristallin.
  - 10. Procédé selon la revendication 9 où le précipitant comprend du polyéthylèneglycol.
  - Procédé selon la revendication 9 où le dimère de M-CSFα est produit à partir de M-CSFα recombiné isolé à partir d'une cellule bactérienne.
  - 12. Procédé selon la revendication 9 où au moins l'un des monomères polypeptidiques est non glycosylé.
  - Procédé selon la revendication 9 où au moins l'un des monomères polypeptidiques contient les résidus 4 à 158 du polypeptide M-CSFα mature.
  - 14. Procédé selon la revendication 9 où le dimère de M-CSF $\alpha$  cristallin est biologiquement actif quand il est resolubilisé.
- 15. Procédé selon la revendication 9 où l'étape de précipitation du cristal de M-CSFα comprend la mise en équilibre du mélange de dimère de M-CSFα avec une solution de précipitant.
  - Procédé selon la revendication 15 où la solution de précipitant comprend une concentration de précipitant plus élevée que le mélange de M-CSFα.
- 55 17. Procédé selon la revendication 15 où l'étape de mise en équilibre comprend l'application du mélange de M-CSFα à une surface et la mise en équilibre du mélange de M-CSFα appliqué avec un réservoir de la solution de précipitant.

- 18. Procédé selon la revendication 9 où l'étape de précipitation du dimère de M-CSFα cristallin comprend une étape d'application de germes cristallins de M-CSFα au mélange de dimère de M-CSFα.
- Procédé selon la revendication 9 où l'étape de précipitation du dimère de M-CSFα comprend l'étape de modification de la température du mélange de M-CSFα.
- 20. M-CSFα (4-158) cristallin.

5

25

30

35

40

45

- 21. M-CSF cristallin selon la revendication 20 comprenant des cristaux de M-CSF ayant des mesures d'au moins environ 0,3 mm dans au moins deux dimensions.
  - 22. M-CSF cristallin selon la revendication 20 qui est capable de diffracter les rayons X pour produire un diagramme de diffraction représentant la structure tridimensionnelle du M-CSF.
- 23. M-CSF cristallin selon la revendication 22 où le diagramme de diffraction représente la structure tridimensionnelle d'une partie du M-CSFα (4-158) avec une précision d'au moins cinq angströms, où la partie inclut les résidus d'acides aminés entre le résidu 4 et le résidu 153 environ.
- 24. Procédé de détermination d'une structure tridimensionnelle d'un dimère de M-CSF contenant deux monomères polypeptidiques de M-CSF ayant chacun entre environ 146 et 162 acides aminés depuis les résidus 1 environ à 5 environ du M-CSF mature, le procédé comprenant:
  - a) la cristallisation du dimère de M-CSF;
  - b) l'irradiation du M-CSF cristallin pour obtenir un diagramme de diffraction caractéristique du M-CSF cristallin; et
  - c) la transformation du diagramme de diffraction en la structure tridimensionnelle du dimère de M-CSF.
  - 25. Procédé d'utilisation d'une structure tridimensionnelle de M-CSF issue d'un cristal de M-CSF où la structure tridimensionnelle du M-CSF inclut une région de liaison de récepteur du M-CSF, le procédé comprenant l'identification de composés ayant des structures qui interagissent avec une région de liaison de récepteur de la structure tridimensionnelle de M-CSF et jouent le rôle d'agoniste ou d'antagoniste du M-CSF.
  - 26. Procédé selon la revendication 25 où la structure tridimensionnelle du M-CSF inclut des coordonnées des carbones alpha sensiblement identiques à celles de l'information structurale présentée dans l'annexe 1.
  - 27. Procédé d'identification d'agonistes ou d'antagonistes du M-CSF, le procédé comprenant les étapes suivantes:
    - a) cristallisation d'un dimère de M-CSF pour former au moins un cristal de M-CSF, le dimère de M-CSF contenant un groupe de résidus d'acides aminés définissant une région de liaison de récepteur du M-CSF;
    - b) irradiation d'un cristal de M-CSF produit par le protocole de l'étape (a) pour obtenir un diagramme de diffraction du cristal de M-CSF;
    - c) détermination de la structure tridimensionnelle du M-CSF d'après le diagramme de diffraction, la structure tridimensionnelle incluant une région de liaison de récepteur du M-CSF; et
    - d) identification d'un composé agoniste ou antagoniste du M-CSF ayant une structure tridimensionnelle qui reproduit fonctionnellement les résidus accessibles aux solvants de liaison de récepteur du M-CSF essentiels présents dans la structure tridimensionnelle de la région de liaison de récepteur du M-CSF, ledit agoniste ou antagoniste du M-CSF ayant une capacité de transduction des signaux modifiée à l'égard des cellules sensibles au M-CSF.
- 28. Procédé selon la revendication 27 où lesdits résidus accessibles aux solvants ne participent pas à la formation de l'interface du dimère de M-CSF.




FIG. 1





FIG. 2



F16.3




FIG. 4A







FIG. 4B










FIG. 6A




FIG. 6B




FIG. 6C

